



MF27/3/13



MF28/3/250/SC



MF29/3/250



C1/12/SC

Bottles



Dreschel bottle head, adjustable

Cat. No.	Socket size	Thread
MF27/3/13	24/29	13

Dreschel bottle heads, plain

Comply with BS 2461.

Cat. No.	Cone size
MF28/3/125/SC	24/29
MF28/3/250/SC	24/29
MF28/3/500/SC	24/29

Dreschel bottle heads, sintered

Comply with BS 2461. Cone 24/29..

Cat. No.	Porosity	Stem length, mm
MF48/03/125	0	130
MF48/03/250	0	165
MF48/03/500	0	215
MF48/13/125	1	130
MF48/13/250	1	165
MF48/13/500	1	215
MF48/23/500	2	215

Dreschel gas washing bottles

Cat. No.	Capacity	Socket size
MF29/3/125	125	24/29
MF29/3/250	250	24/29
MF29/3/500	500	24/29

Condensers



With screw thread connections (except air condensers). All condensers are supplied complete with screw thread connectors (4510/02).

Liebig

A condenser designed to give turbulent flow in the water jacket with a comparatively slow rate of flow of coolant water by use of a small annular space.

Cat. No.	Approx surface area x 10 ² m ²	Socket size	Cone size	Approx. effective x overall length, mm
C1/11/SC	0.60	14/23	14/23	155 x 260
C1/12/SC	1.00	19/26	19/26	208 x 340
C1/13/SC	1.00	24/29	24/29	250 x 390
C1/22/SC	2.00	19/26	19/26	395 x 530
C1/23/SC	2.00	24/29	24/29	395 x 535
Spares				
FF/16202	Spare screw thread tubing connector			
FF/16702	Spare screw thread tubing connector PTFE			