



EX5/55/100

Extraction liquid-solid, Soxhlet extractors

To BS 2071 where applicable. Designed primarily for use with proprietary makes of paper thimble and suitable sizes are listed on page 203. Large samples may be placed directly in the extractor body.

Cat. No.	Approx surface area x 10 ² m ²	Socket size	Cone size	Approx. effective x overall length, mm
EX5/21	6	19/26	14/23	150
EX5/32	20	24/29	19/26	175
EX5/43	40	29/32	24/29	200
EX5/53	60	34/35	24/29	205
EX5/55	60	34/35	34/35	205
EX5/55/100	100	34/35	34/35	240
EX5/63	100	40/38	24/29	220
EX5/83	200	50/42	24/29	250
EX5/75	600	75mm*	34/35	370
EX5/105	2000	100mm*	34/35	535

* Approximate height is measured from the wide end of the cone.



JC75F

Adapters and clips

Items EX5/75 and EX5/105 are designed to accommodate large samples. They utilise 75mm and 100mm at flange joints and require the at socket adapters and clips listed below.

Cat. No.	Description	To fit	Socket size	Flange size, mm
DA5/75	Adapter	EX5/75	34/35	75
DA5/100	Adapter	EX5/105	34/35	100
JC75F	Clip	EX5/75	-	75
JC100F	Clip	EX5/105	-	100



JC100F

Complete standard assemblies.

Cat. No.	Nominal Capacity, ml	Extractor etc.	Components Flask	Condenser
EX/12008	100	EX5/63	FR250/3S	CX7/06/SC
EX/12010	200	EX5/83	FR500/3S	CX7/08/SC
EX/12012	600	EX5/75 DA5/75 JC75F	FR2L/5S	CX5/25/SC



EX/12012