

Flasks, round bottom, short neck

Cat. No.	Nominal cap, ml	Socket size	Approx overall height, mm	Approx. neck length, mm
FR5/0S	5	10/19	50	5
FR5/1S	5	14/23	65	15
FR10/0S	10	10/19	65	15
FR10/1S	10	14/23	70	15
FR25/1S	25	14/23	75	10
FR25/2S	25	19/26	80	15
FR25/3S	25	24/29	80	15
FR50/1S	50	14/23	90	15
FR50/2S	50	19/26	90	15
FR50/3S	50	24/29	95	15
FR50/5S	50	34/35	95	10
FR100/1S	100	14/23	105	15
FR100/2S	100	19/26	108	15
FR100/3S	100	24/29	110	15
FR100/4S	100	29/32	114	15
FR100/5S	100	34/35	112	10
FR150/2S	150	19/26	114	15
FR150/3S	150	24/29	115	15
FR150/4S	150	29/32	120	15
FR150/5S	150	34/35	118	10
FR250/1S	250	14/23	126	15
FR250/2S	250	19/26	129	15
FR250/3S	250	24/29	130	15
FR250/4S	250	29/32	135	15
FR250/5S	250	34/35	138	15
FR500/2S	500	19/26	147	15
FR500/3S	500	24/29	150	15
FR500/5S	500	34/35	166	25
FR1L/2S	1000	19/26	177	15
FR1L/3S	1000	24/29	180	15
FR1L/5S	1000	34/35	196	25
FR1L/7S	1000	45/40	201	25
FR2L/3US	2000	24/29	205	10
FR2L/3S	2000	24/29	215	20
FR2L/4S	2000	29/32	224	20
FR2L/5S	2000	34/35	228	25
FR3L/3S	3000	24/29	240	20
FR3L/5S	3000	34/35	249	25
FR4L/7S	4000	45/40	290	25
FR5L/4S	5000	29/32	270	20
FR5L/5US	5000	34/35	275	20
FR6L/7S	6000	45/40	300	25
FR10L/5S	10000	34/35	350	25



Short neck

Flasks, round bottom, medium neck

Cat. No.	Nominal cap, ml	Socket size	Approx overall height, mm	Approx. neck length, mm
FR50/3M	50	24/29	130	45
FR100/3M	100	24/29	145	45
FR150/3M	150	24/29	150	45
FR250/3M	250	24/29	165	45
FR500/3M	500	24/29	180	45



Medium neck