



SB stoppers

## Stoppers

Hollow blown to reduce weight.

Cat. No.	Cone size
SB10	10/19
SB12	12/21
SB14	14/23
SB19	19/26
SB24	24/29
SB29	29/32
SB34	34/35



SH4A

## Fixed immersion thermometer pocket

Closed end, made to fit Quickfit stillheads. It will take a normal unjointed laboratory thermometer.

Cat. No.	Cone size	Approx. stem length, mm
SH4A	14/23	45

## Test tubes

With Quickfit® precision ground joint. Accepts stoppers on page 405.  
All dimensions and capacities are approximate.



MF24/3

Cat. No.	Socket size	Overall, mm	Dia., mm	Cap., ml
MF24/0/4	10/19	100	14	5
MF24/12/4	12/21	100	17	10
MF24/1	14/23	100	19	11
MF24/1/5	14/23	125	19	15
MF24/1/6	14/23	150	19	19
MF24/2	19/26	100	23	19
MF24/2/6	19/26	150	23	31
MF24/3	24/29	150	29	54
MF24/3/8	24/29	200	29	75
MF24/5	34/35	250	41	195




BC24

## Centrifuge tubes

Each tube is marked with its weight to the nearest 0.1g to enable tubes of similar weight to be selected from any batch. Max RCF 3000.

Cat. No.	Socket size	Overall, mm	Dia., mm	Cap., ml
BC24/C14T	14/15	17 to 17.3	115	10
BC24/C24R	24/20	28 to 29	115	30



MF23

## Drying, angled

Cat. No.	Cone size
MF23/1	14/23
MF23/2	19/26
MF23/3	24/29
MF23/4	29/32