



TR/36

Test tube peg racks, filled PP

NALGENE®

Suitable for sample collection tubes, reaction tubes or similar-diameter tubes. Will not float in a water bath. Can be used to dry inverted tubes. Also can be used to hold and dry electrophoresis and chromatography plates. Chemical resistant. (Autoclavable)

Cat. No.	Code	Tube sizes	(lxwxh), mm	Array	Pk Qty
White racks					
TR/36040	5977-0013	10-13 mm	179x127x64	8x12	2
TR/36042	5977-0017	14-17 mm	187x105x70	5x10	2
Blue racks					
TR/36340	5977-0313	10-13 mm	179x127x64	8x12	2
TR/36342	5977-0317	14-17 mm	187x105x70	5x10	2

Test tube racks, Unwire™, ResMer™

NALGENE®

Unwire test tube racks securely hold a full load of test tubes or centrifuge tubes in the following popular sizes: 13, 16, 20, 25 and 30mm. Up to 72 of the smaller diameter tubes can be held at once. All sizes are moulded in six bright, permanent colours. ResMer™ manufacturing technology makes these racks keep their bright colours longer. Useful in a wide temperature range with good chemical resistance. (Autoclavable).



TR/36 Unwire

Cat. No.	Code	Tube sizes	(lxwxh), mm	Array
White racks				
TR/36002	5970-0013	13	200x102x57	6x12
TR/36004	5970-0016	16	248x127x70	6x12
TR/36006	5970-0020	20	250x102x83	4x10
TR/36008	5970-0025	25	300x121x92	4x10
TR/36010	5970-0030	30	283x108x83	3x8
Orange racks				
TR/36102	5970-0113	13	200x102x57	6x12
TR/36104	5970-0116	16	248x127x70	6x12
TR/36106	5970-0120	20	250x102x83	4x10
TR/36108	5970-0125	25	300x121x92	4x10
TR/36110	5970-0130	30	283x108x83	3x8
Yellow racks				
TR/36202	5970-0213	13	200x102x57	6x12
TR/36204	5970-0216	16	248x127x70	6x12
TR/36206	5970-0220	20	250x102 x83	4x10
TR/36208	5970-0225	25	300x121x92	4x10
TR/36210	5970-0230	30	283x108x83	3x8
Blue racks				
TR/36302	5970-0313	13	200x102x57	6x12
TR/36304	5970-0316	16	248x127x70	6x12
TR/36306	5970-0320	20	250x102x83	4x10
TR/36308	5970-0325	25	300x121x92	4x10
TR/36310	5970-0330	30	283x108x83	3x8
Green racks				
TR/36402	5970-0413	13	200x102x57	6x12
TR/36404	5970-0416	16	248x127x70	6x12
TR/36406	5970-0420	20	250x102x83	4x10
TR/36408	5970-0425	25	300x121x92	4x10
TR/36410	5970-0430	30	283x108x83	3x8
Red racks				
TR/36502	5970-0513	13	200x102x57	6x12
TR/36504	5970-0516	16	248x127x70	6x12
TR/36506	5970-0520	20	250x102x83	4x10
TR/36508	5970-0525	25	300x121x92	4x10
TR/36510	5970-0530	30	283x108x83	3x8