

More space-efficient than full-size test tube racks. Fit easily in ice or water baths and won't float. All the features of full-size Unwire test tube racks. ResMer™ manufacturing technology makes these racks keep their bright colours longer. Useful in a wide temperature range with good chemical resistance. (Autoclavable).

Cat. No.	Code	Tube sizes	(lxwxd), mm	Array
White racks				
TR/36012	5972-0013	13 mm	102x102x56	6x6
TR/36014	5972-0016	16 mm	126x126x68	6x6
TR/36016	5972-0020	20 mm	128x103x83	4x5
TR/36018	5972-0025	25 mm	122x122x75	4x4
TR/36020	5972-0030	30 mm	109x109x84	3x3
Blue racks				
TR/36312	5972-0313	13 mm	102x102x56	6x6
TR/36314	5972-0316	16 mm	126x126x68	6x6
TR/36316	5972-0320	20 mm	128x103x83	4x5
TR/36318	5972-0325	25 mm	122x122x75	4x4
TR/36320	5972-0330	30 mm	109x109x84	3x3
Green racks				
TR/36412	5972-0413	13 mm	102x102x56	6x6
TR/36414	5972-0416	16 mm	126x126x68	6x6
TR/36416	5972-0420	20 mm	128x103x83	4x5
TR/36420	5972-0430	25 mm	109x109x84	3x3
Red racks				
TR/36512	5972-0513	13 mm	102x102x56	6x6
TR/36514	5972-0516	16 mm	126x126x68	6x6
TR/36516	5972-0520	20 mm	128x103x83	4x5
TR/36520	5972-0530	25 mm	109x109x84	3x3



TR/36 Unwire

racks, test tube

Test tube stands, wooden, plain

Hardwood with plastic end supports. Packed flat for transport and storage.

Cat. No.	Number of holes x size
TR/47002	6 x 22 mm
TR/47004	12 x 22 mm
TR/47006	6 x 32 mm
TR/47008	2 x 32 mm & 4 x 22 mm



TR/4700 and TR/4701

Test tube stands, wooden, with drying pegs

Hardwood with plastic end supports.

Cat. No.	Number of holes x size
TR/47012	6 x 22 mm
TR/47014	12 x 22 mm
TR/47018	2 x 32 mm & 4 x 22 mm

Test tube racks

Clifton Range®

Test tube racks suitable for use with Clifton baths. Made of high grade stainless steel. Overall dimensions of rack: 270mm long x 70mm wide, height to top of lifting handle: 138mm. Max number of racks per bath:-
Bath size 4 litre=1 rack, 8 litre=2 racks, 14 litre=4 racks, 22/28 litre=6 racks

Cat. No.	Code	Description
BW/85250	6872	Test tube rack, stainless steel, 36 holes x 13mm dia
BW/85252	6870	Test tube rack, stainless steel, 26 holes x 17mm dia
BW/85254	6873	Test tube rack, stainless steel, 18 holes x 19mm dia
BW/85256	6871	Test tube rack, stainless steel, 16 holes x 26mm dia
BW/85258	6900	Test tube rack, stainless steel, 12 holes x 32mm dia



BW/85250