




TR/84002



TR/930



ST/13002

Test tube rack, Z pattern, aluminium

Z pattern test tube rack No of Holes x Diameter (mm) 8 x 13mm, 6 x 19mm and 5 x 28mm. Height 76mm x Length 203mm x Width 114mm.

Cat. No.	Description
TR/84002	Z pattern test tube rack

Test tube racks, nylon coated mild steel

Cat. No.	Tube Size	Format
TR/93002	13mm	4 x 6
TR/93004	13mm	3 x 12
TR/93006	16mm	3 x 4
TR/93008	16mm	4 x 6
TR/93010	16mm	3 x 12
TR/93012	20mm	3 x 6
TR/93014	20mm	3 x 12
TR/93016	25mm	3 x 4
TR/93018	Universal bottles	2 x 6

Stands, petri dish

- Holds up to six 90mm petri dishes
- Manufactured from epoxy-coated steel
- Autoclavable
- Ideal for small batch incubation and moving samples

Azlon

Cat. No.	Code	Colour
ST/13002	SWP062	White

Refractometers

Refractometer, Kern

KERN

A range of economically priced refractometers in 3 versions: hand-held, digital and bench type.

The analogue hand-held version is universal, maintenance free; the optical system and prism cover are manufactured from a special material that allows low-tolerance measuring.

The digital version is again maintenance free, has a large TFT colour digital display with integrated temperature reading. Manufactured to international IP65 protection class against dust, and water splashes.

The Abbe bench type refractometer allows easy, efficient and sustainable use in everyday applications. An integrated scale allows the unit to be used for different applications and permits reading of measurement results accurately.

Please note: Other models are available for different applications and ranges, please enquire.



RF/25100



RF/25120

Technical Specifications

Model	ORA 1RE	ORF 1RS	ORT 1RS
Scale	Refractive index	Refractive index	Refractive index/Brix
Measuring range	1,333 - 1,405 nD 1,405 - 1,468 nD 1,468 - 1,517 nD	1,3330 - 1,5400 nD	Brix 0-95% R.I. 1,3000 - 1,7000 nD
Division	0,005 nD	±0,0001 nD	Brix 0,25% R.I. 0,0005 nD
Accuracy	---	±0,0005 nD	Brix ± 0,1% R.I. ± 0,0002 nD

Cat. No.	Code	Description
RF/25100	ORA 1RE	Refractometer, hand-held 1,333 - 1,517 nD
RF/25120	ORF 1RS	Refractometer, digital RI 1,3330 - 1,5400 nD
RF/25140	ORT 1RS	Refractometer, bench type, Abbe RI 1,3000 - 1,7000 nD, Brix 0 - 95%