



RF/410

Certified Reference Materials for Refractometers



- Aqueous solutions up to 40 °Brix
- Calibration Oils for high range applications
- High accuracy
- Traceable to NIST and ICUMSA
- UKAS certified
- 12-months validity

AG Fluids - Multi-pack of 5 x 5ml bottles

Cat. No.	Product	Code	Refractive Index	°Brix
RF/41012	AG2.5	90-501	1.33659	2.50
RF/41014	AG5	90-502	1.34026	5.00
RF/41016	AG7.5	90-503	1.34401	7.50
RF/41018	AG10	90-504	1.34782	10.00
RF/41020	AG11.2	90-505	1.34968	11.20
RF/41022	AG12	90-506	1.35093	12.00
RF/41024	AG12.5	90-507	1.35171	12.50
RF/41026	AG15	90-508	1.35568	15.00
RF/41028	AG40	90-518	1.39986	40.00

Maximum Uncertainty: ±0.000037 RI ±0.019 °Brix

Calibration Oil - Multi-pack of 5 x 5ml bottles

Cat. No.	Code	Refractive Index	°Brix
RF/41050	90-525	1.46990	71.81

Maximum Uncertainty: ±0.000074 RI; ±0.030 °Brix



RF/75802

Refrigerators

Refrigerators and freezers, basic



Labcold basic fridges and freezers are available in two sizes only: large upright and small under counter. Manufactured to the high standard you would expect from Labcold, these fridges and freezers are laboratory standard at a budget price, featuring robust enameled exteriors and wipe clean plastic interiors.

Cat. No.	Model	Description	Ext Dimensions (HxWxD) mm	Capacity litres	Shelves
RF/75802	RLPL04043	Refrigerator, basic, under counter 105 litres	841x495x500	105	3+1 drawer
RF/75804	RLFR1804	Refrigerator, basic, upright 505 litres	1710x765x725	505	4
RF/75806	RLVL03203	Freezer, basic, under counter 70 litres	835x500x490	70	2+3 drawers
RF/75808	RLVF1825	Freezer, basic, upright 505 litres	1710x765x725	505	6

Notes
