

Refrigerators, Spark Free



Labcold Spark-free refrigerators are manufactured with an interior having all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it a significant aid to laboratory safety. For extra security a built in lock is fitted.

The temperature display features high and low audio and visual temperature alarms along with a door open alarm. The controller also records the minimum and maximum temperature since the last reset. For improved temperature stability a fan is included for more even circulation and air distribution. As an added feature, off cycle automatic defrosting is included, complete with 6 adjustable shelves.

NB: Not suitable for the storage of open containers of volatile substances

Technical Specifications

Capacity (litres)	400
Temperature range	0°C to +10°C
Factory set point	+4°C
Exterior dimensions (HxWxD)	1855 x 600 x 585mm
Chamber dimensions (HxWxD)	1620 x 506 x 485mm
Weight	69Kg
Shelves	6
Shelf dimensions (WxD)	5 x 500 x 410mm plus 1 x 500 x 230mm
Door Type	Solid
Lock Type	Key lock with two keys

An access port is available as an optional extra, please enquire.

Cat. No.	Model	Description
RF/75022	RLPR1514	Spark-free refrigerator, 400 litres



RF/75022

Labcold specialises in supplying the medical and scientific community with refrigerators and freezers designed to meet the requirements of the most demanding laboratory.

Our range includes:

- Sparkfree fridges & freezers
- Blood product storage
- Pharmacy refrigeration
- Reagent fridges
- ULT freezers
- Cooled incubators

