



RF/75308

Freezer, Spark-free, upright



Labcold Spark-free upright freezer, manufactured with an interior having all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it a significant aid to laboratory safety. For extra security a built in lock is fitted. The temperature display features high and low audio and visual temperature alarms along with a door open alarm. The controller also records the minimum and maximum temperature since the last reset, defrost is actioned manually, supplied with 6 fixed shelves.

NB: Not suitable for the storage of open containers of volatile substances

Technical Specifications	
Capacity (litres)	400
Temperature range	-12°C to -24°C
Factory set point	-20°C
Exterior dimensions (HxWxD)	1855 x 600 x 585mm
Chamber dimensions (HxWxD)	1620 x 470 x 480mm
Weight	79Kg
Shelves	6
Shelf dimensions (WxD)	480 x 410mm
Door Type	Solid
Lock Type	Built in

An access port is available as an optional extra, please enquire.

Cat. No.	Model	Description
RF/75308	RLVF1514	Spark-free freezer, 400 litres



RF/75508

Freezers, chest



Completely spark-free, the chest configuration of these freezers helps reduce temperature rises when the door is opened.

- Spark-free interior
- Door lock
- Easy-grip door handle
- CFC, HCFC and ammonia free
- Temperature set-point range: -18 to -25°C
- Visual high temperature alarm
- Operator controlled defrost

Cat. No.	Model	Ext. dimensions hwxwd (mm)	Capacity litres	Storage baskets
RF/75502	RLCF0720	860x750x730	215	1
RF/75508	RLCF1520	860x1300x730	447	2
RF/75510	RLCF2120	860x1700x730	607	3

Notes