

Safety clothing

Apron, polythene, disposable – flat dispenser pack

Disposable polythene aprons, with halter neck and waist ties are supplied in convenient, at wall mountable or worktop dispenser packs which allow easy, non-stick removal of each apron. Resistant to fluid and bacteria penetration.

| Cat. No. | Size (w x l) cm | Colour | Inner | Case |
|----------|-----------------|--------|-------|------|
| SC/66002 | 69 x 107 | White | 100 | 1000 |
| SC/66004 | 69 x 122 | White | 100 | 1000 |

Apron, polythene, disposable - on a roll

Aprons on a roll are manufactured from the same material as the dispenser pack version but are supplied on rolls for convenience. 69 x 117cm - Longer length available

| Cat. No. | Colour | Inner Pk | Case |
|-----------|----------------------------|----------|------|
| SC/66102 | White | 200 Roll | 1000 |
| SC/66104 | Blue | 200 Roll | 1000 |
| SC/66106 | Red | 200 Roll | 1000 |
| Accessory | | | |
| SC/66114 | Dispenser, plastic covered | Each | Each |

Aprons, red rubber

High quality, flexible and waterproof rubber apron for heavy duty applications. Hard wearing, offering excellent resistance to fats, oils and high temperatures. One size to fit all: 91 x 107cm (36 x 42"). Supplied with adjustable ties. Each apron individually wrapped in a polybag, 20 pieces per carton.

| Cat. No. | Description |
|----------|-------------------------------|
| SC/66120 | Apron, red rubber, 91 x 107cm |

Apron, cryo

Cryo aprons, blue, with multi-layers to protect the torso and legs against temperature from -160°C to + 150°C. The waterproof nylon outer cover offers protection from contact with ultra- cold surfaces and from splashes. Ideal for working in deep-freeze environments, with cryogenic liquids, dry ice, liquid nitrogen (vapour phase) or when working with autoclaves. The aprons have adjustable straps and quick release buckles at the waist and neck.

| Cat. No. | Description |
|----------|-------------------------------------|
| CR/15000 | Apron, cryogenic, blue, 915mm long |
| CR/15002 | Apron, cryogenic, blue, 1070mm long |
| CR/15004 | Apron, cryogenic, blue, 1220mm long |
| CR/15006 | Apron, cryogenic, blue, 1380mm long |

Laboratory coats

Visitor's lab-coat, white, breathable, lightweight material. Designed for light duty applications and offers a barrier against non-hazardous particles and liquid splashes, with a zip fastener, available in sizes small, medium, large, extra-large, 2XL, 3XL and 4XL*.

| Cat. No. | Description |
|----------|--|
| SC/68320 | Visitor's lab coat, lightweight breathable |

Coverall, white, designed to offer cost effective protection against hazardous dusts and limited Type 5/6 liquid splashes. With elasticated hood, ankles, cuffs and waist, 2 way zip with storm flap, ultra-low linting, Available in sizes small, medium, large, extra-large, 2XL, 3XL and 4XL*.

| Cat. No. | Description |
|----------|--|
| SC/68324 | Lab coverall, protects against dust and limited Type 5/6 liquid splashes |

*Please state size when ordering



SC/661



SC/66120



CR/1500



SC/68320

SC/68324