



SG/66008

Gloves, heat protection, 'Crusader Flex'

Nitrile dip coating on a felt liner. Heat protective coating that extends the entire length of the glove, both front and back. Hard wearing. Designed for intermittent handling of medium to hot objects up to 200°C. Excellent grip. (Sold as a pair). Length, 330mm

Cat. No.	Code	Description	Size
SG/66008	42-474	Crusader Flex	10



SG/28204

Thermal Glove, Navy

Can be used for protection to -25°C.

Cat. No.	Size
SG/28204	One Size



CR/9110

Gloves, Cryogenic

Designed to provide protection down to -190°C. Supplied in pairs.

Cat. No.	Description
CR/91100	Cryo gloves 7½ - 8
CR/91102	Cryo gloves 8½ - 9
CR/91104	Cryo gloves 9½ - 10



SH/1600

Hand/finger protectors, silicone rubber

- Ideal for the safe handling of hot beakers, flasks, etc.
- Capable of withstanding temperatures up to 250°C.
- Also capable of handling cold objects down to -57°C.
- Contact surfaces are multi-studded for positive grip.

Cat. No.	Reference	Description
SH/16002	PGR402	Medium
SH/16006	PGR406	Finger