

Safety storage cabinets

Systemlabor FWF90 safety cabinets are designed to ensure the safe storage of flammable chemicals. The cabinets are re rated to 90 minutes at 200°C and certified in accordance with DIN 12925-1. In the event of fire they keep the temperature of the stored chemicals below critical levels, long enough to allow safe evacuation of the laboratory. Acid storage cabinets are constructed of solid propylene for the storage of corrosive chemicals. All safety cabinets are ventilated to avoid the build-up of noxious fumes. Please ask for details of the full range of products and a quotation.

Hazardous substances and flammable liquid storage cabinets

- Flush fitting handles
- Adjustable perforated shelves
- Liquid tight sumps
- Robust steel construction
- Reinforced doors
- Easy drain perforated shelves (galvanised)
- Removable sump tray
- Hi-Viz warning yellow colour
- Warning symbols supplied
- Tough powder coated finish
- Lockable

Cat. No.	Code	Description	Overall dimensions (HxWxD) mm	Shelves
SC/24262	HZ5	Single door cabinet	700X350X300	1
SC/24264	HZ6	Single door cabinet	700X900X460	1
SC/24266	HZ3	Single door cabinet	900X460X460	1
SC/24268	HZ4	Double door cabinet	900X900X460	1
SC/24270	HZ1	Double door cabinet	1800X900X460	3



SC/242



SC/24262

Acid and alkali storage cabinets

These cabinets enable safe and identifiable storage in the workplace

- Strong steel construction
- Easy-drain perforated shelves (galvanised)
- Removable sump
- Durable clean white finish
- Resilient powder coating
- Lockable

Cat. No.	Code	Description	Overall dimensions (HxWxD) mm	Shelves
SC/24142	AA5	Acid & alkali storage cabinet	700x350x300	1
SC/24144	AA6	Acid & alkali storage cabinet	700X900X460	1
SC/24146	AA3	Acid & alkali storage cabinet	900x460x460	1
SC/24148	AA4	Acid & alkali storage cabinet	900x900x460	1
SC/24150	AA1	Acid & alkali storage cabinet	1800x900x460	3



SC/241



SC/24142