



SC/24304

Fireproof safety cabinets – EN14470-1 Type 30min – FM approved

Fireproof safety cabinet, 30 min fire-rated for storing flammable substances! Designed, tested and approved according to EN14470-1 standard that requires the safest features to protect your flammable chemical from fire. Also available in 90 and 60 minutes.

- Tested and approved according to EN14471-1 and EN1363-1 and FM
- 30 minutes fire resistant built in HPL
- Self-closing door(s) with thermo-fuse
- Removable retention tank
- Adjustable retention shelves
- Normalized label
- Ventilation inlet and outlet with thermo-fuse
- Pedestal area to move easily with a forklift
- Doors with cylinder locks

Cabinets to store up to 4 different substances are also available.

Technical Specifications

3031E	W605 x D630mm x H635mm – supplied with 1 base trough and 1 perforated cover plate
3032E	W1115mm x D630mm x H635mm – supplied with 1 base trough and 1 perforated cover plate
3034E	W605mm x D550mm x H1980mm – supplied with 1 base trough and 3 shelves
3035E	W1115mm x D550mm x H1980mm – supplied with 1 base trough and 3 shelves
3032M11	W1115mm x D630mm x H635mm – 2 compartments – supplied with 2 base troughs and 2 perforated cover plates
3034M2	W605mm x D550mm x H1980mm – 2 compartments - supplied with 2 base troughs and 4 shelves
3035M11	W1115mm x D550mm x H1980mm – 2 compartments - supplied with 2 base troughs and 6 shelves
3035M22	W1115mm x D550mm x H1980mm – 4 compartments - supplied with 4 base troughs and 8 shelves

Cat. No.	Model	Description
SC/24300	3031E	Ecosafe under bench cabinet – one door
SC/24302	3032E	Ecosafe under bench cabinet – two doors
SC/24304	3034E	Ecosafe tall cabinet – one door
SC/24306	3035E	Ecosafe tall cabinet – two doors
SC/24308	3032M11	Ecosafe multi-risk under bench cabinet – two doors – two compartments
SC/24310	3034M2	Ecosafe multi-risk tall cabinet – one door – two compartments
SC/24312	3035M11	Ecosafe multi-risk tall cabinet – two doors – two compartments
SC/24314	3035M22	Ecosafe multi-risk tall cabinet – two doors – four compartments
Accessories		
SC/24316	E48	Retention shelf for one door cabinets
SC/24318	E35	Retention shelf for two doors cabinets
SC/24320	T148	Drawer for one door cabinets
SC/24322	T235	Drawer for two doors cabinets
SC/24324	C148	Perforated cover plate for one door cabinet
SC/24326	C235	Perforated cover plate for two door cabinet
SC/24328	CDF-ACORG	Filtration Box
SC/24330	CDV-A	Ventilation Box
SC/24332	KRC	Connection Kit



SC/24304