

Safety, waste disposal

Waste bin, pedal operated

White plastic mobile pedal bin, capacity 60L. Overall dimensions - 700 x 380 x 490mm.

Cat. No.	Description
SD/91002	Pedal bin



SD/91002

Disposal bins, cardboard & LDPE

- Self assembly rigid boxes for safe disposal of broken glass.
- Glass is safely held in strong, leak-proof and puncture resistant LDPE liner.
- Lid has an integral waste flap for easy use and has a sealing feature for use when full.
- Boxes are clearly printed with instructions for use.

Azlon

Cat. No.	Code	Description	Size (WxDxH, cm)
SD/13002	DBG001P	Glass disposal bin floor standing	30 x 30 x 70
SD/13004	DBG002P	Glass disposal bin bench top	23 x 23 x 26.5



SD/1300

Sharps containers – Sharpsafe

The Sharpsafe containers are available in various sizes and comply with BS 7320:1990 and UN3291. Also available in the range is the Clinisafe, a multi-purpose clinical waste collection unit for Cytotoxic Placental and Laboratory waste and the Quiver, an ideal container for holding laparoscopic instruments.

Cat. No.	Description	Pk Qty
SD/26402	Sharpsafe 2 litre	50
SD/26404	Sharpsafe 4 litre	50
SD/26406	Sharpsafe 7 litre	50
SD/26408	Sharpsafe 13 litre	20
SD/26410	Sharpsafe 24 litre	10
SD/26412	Sharpsafe pocket size	100
SD/26416	Quiver 7.5 litre	20
SD/26418	Clinisafe 9 litre	20



SD/2640

Disposal bags

Biohazard autoclave bags with warning instruction and sterilization indicator field, strong 50µm Polypropylene film with broad seam at bottom. Autoclavable at 134°C (When autoclaving, the bag must not be completely sealed).

Cat. No.	Dimensions (WxDxH, mm)	Volume (approx) Litres	Pk Qty
AC/16300	300 x 500	8	500
AC/16302	400 x 780	38	500
AC/16304	600 x 800	60	500
AC/16306	700 x 1100	110	350



AC/1630