



IS/162



BS/1220



SC/8400



SF/9300



SF/93010

Scoops

Scoops, measuring, polypropylene



- Made from opaque white polypropylene
- Top is flat for levelling of powders
- Base is flat enabling scoops to stand horizontally
- The capacity is clearly moulded on each size

Cat. No.	Code	Capacity	Overall length	Pk Qty
IS/16218	SWL006	5ml	82mm	5
IS/16228	SWL008	10ml	102mm	5
IS/16234	SWL012	25ml	137mm	5
IS/16240	SWL018	50ml	160mm	5
IS/16246	SWL024	100ml	198mm	5
IS/16256	SWL030	250ml	259mm	5
IS/16266	SWL036	500ml	315mm	5
IS/16276	SWL140	1000ml	385mm	1

Scoops, weighing, borosilicate glass



Ideal for weighing small quantities of powder and then introducing into volumetric flasks etc.

Cat. No.	Code	Capacity	Length	Pk Qty
BS/12202	2115/02M	3ml	70mm	10
BS/12204	2115/04M	6ml	85mm	10
BS/12206	2115/06M	10ml	100mm	10

Scoops, cylindrical, stainless steel

Cat. No.	Capacity	Diameter x overall length
SC/84002	100ml	80 x 180mm
SC/84004	250ml	100 x 260mm
SC/84006	500ml	100 x 325mm
SC/84008	1000ml	115 x 340mm

Sealing film

Sealing film, Parafilm M

Designed to stretch over and seal a wide range of laboratory vessels. Supplied in dispensing boxes containing rolls with paper backing to prevent sticking. Parafilm is not suitable for use with solvents.

Cat. No.	Roll width x length
SF/93002	50mm x 75m (2 x 246ft)
SF/93004	100mm x 38m (4 x 124ft)
SF/93006	100mm x 75m (4 x 248ft)
SF/93008	500mm x 15m (19.6 x 49ft)

Parafilm cutter

Sturdy transparent acrylic dispenser. Protected cutter with very sharp blade. Holds 1 or 2 rolls of the 50mm width film or 1 of the 100mm width.

Cat. No.	Description
SF/93010	Parafilm cutter