

Baths, advanced, SUB Aqua Pro series

Grant

The SUB Aqua Pro range is built to the highest standards and specifications and incorporates the latest technology to provide accurate temperature control in even the most demanding applications in research, pharma and biotech. Choose from eight models, including shallow and dual baths with stainless steel tanks and digital controls.

- Ambient +5 to 99°C
- Stability $\pm 0.2^\circ\text{C}$
- Set and Forget™ technology for fast heat up and accurate temperature control
- Advanced dry start and run dry protection
- Comprehensive safety features to protect your samples and your workspace
- 3-year warranty as standard
- 2-point user calibration
- 3 programmable temperature pre-sets
- Countdown timer (1-999 minutes)
- Polycarbonate lid and base tray included



BW/65646 with lid

Technical Specifications

Model	SAP2	SAP2S	SAP5	SAP12
Tank capacity	2L	2L shallow	5L	12L
Temperature range	ambient +5 to 99°C			
Temperature setting resolution	0.1°C			
Stability @ 70°C	$\pm 0.2^\circ\text{C}$			
Temperature setting/energy regulation	digital			
Temperature display	3 digit bright LED			
Working area (lxwx d mm)	131 x 117 x 132	289 x 139 x 47	281 x 131 x 132	306 x 281 x 132
Drain tap included	No	No	No	Yes
Power consumption 230V/120V/kW	0.25/0.25	0.35/0.35	0.35/0.35	0.8/0.8
Supply voltage	120V or 230V			
Sample protection	Adjustable over temperature cut-off and alarm			

Technical Specifications

Model	SAP18	SAP26	SAP34	SAPD
Tank capacity	18L	26L	34L	5L and 12L
Temperature range	ambient +5 to 99°C			
Temperature setting resolution	0.1°C			
Stability @ 70°C	$\pm 0.2^\circ\text{C}$			
Temperature setting/energy regulation	digital			
Temperature display	3 digit bright LED			
Working area (lxwx d mm)	281 x 485 x 132	278 x 481 x 182	281 x 635 x 180	281 x 131 x 132 and 306 x 281 x 132
Drain tap included	Yes	Yes	Yes	Yes
Power consumption 230V/120V/kW	1.4/1.05	1.4/1.05	1.8/1.3	1.15/1.15
Supply voltage	120V or 230V			
Sample protection	Adjustable over temperature cut-off and alarm			

Cat. No.	Code	Description
BW/65640	SAP2	Water bath, SUB Aqua Pro 2L
BW/65642	SAP2S	Water bath, SUB Aqua Pro 2L shallow
BW/65644	SAP5	Water bath, SUB Aqua Pro 5L
BW/65646	SAP12	Water bath, SUB Aqua Pro 12L
BW/65648	SAP18	Water bath, SUB Aqua Pro 18L
BW/65650	SAP26	Water bath, SUB Aqua Pro 26L
BW/65652	SAP34*	Water bath, SUB Aqua Pro 34L
BW/65654	SAPD	Water bath, SUB Aqua Pro Dual 5L and 12L

*Supplied with a stainless steel lid as standard, not Polycarbonate



BW/65654