



ST/93002

## Stands

### Stand, pipette, polythene covered wire

Horizontal stand for up to 10 pipettes. Dimensions 125 x 156 x 118mm.

| Cat. No. | Description   |
|----------|---------------|
| ST/93002 | Pipette stand |



PS/12034

### Stand, pipette, vertical

White chemically resistant polypropylene stand. Very stable with 230mm diameter base plate. Holds up to 94 pipettes vertically.

| Cat. No. | Code   | Description                       |
|----------|--------|-----------------------------------|
| PS/12034 | SWP033 | Pipette stand, vertical, holds 94 |



ST/9410

### Retort bases, steel

With 10mm tapped bush in centre of short side. All bases have rubber feet. (Supplied without rods)

| Cat. No.                       | Dimensions, mm |
|--------------------------------|----------------|
| Finished in blue acrylic paint |                |
| ST/94102                       | 160 x 100      |
| ST/94104                       | 200 x 125      |
| ST/94106                       | 250 x 160      |
| ST/94108                       | 315 x 200      |



ST/9430

### Rods

All 10mm threaded.

| Cat. No.   | Cat. No.  | Cat. No.        | Dimensions, diam. x length |
|------------|-----------|-----------------|----------------------------|
| Mild steel | Aluminium | Stainless steel |                            |
| ST/94302   | ST/94402  | ST/94522        | 10 X 500mm                 |
| ST/94304   | ST/94404  | ST/94524        | 12.5 X 600mm               |
| ST/94306   | ST/94406  | ST/94526        | 12.5 X 750mm               |
| ST/94308   | ST/94408  | ST/94528        | 12.5 X 1000mm              |



ST/44216 &amp; ST/94008

### Bossheds, alloy

Suitable for rods up to 16mm diameter.

| Cat. No. | Description  |
|----------|--|
| ST/44216 | With offset jaws and nickel plated brass tommy screws. |
| ST/94008 | With blue plastic headed nickel plated steel screws.   |