


Clamps, aluminium alloy

Will accept articles from 2 – 90mm.

Cat. No.	Description
ST/94010	Clamp with cork lined jaws
ST/94012	3-prong clamp with rubber covered jaws



ST/9401

Ring with boss-head, eccentric series

Constructed in Zinc plated mild steel with welded joints, fitted with an alloy boss-head which is die-cast on to the stem and fitted with blue plastic headed nickel steel screws.

These rings have a standard stem length which measures approximately 60mm to the centre of the boss-head.

Cat. No.	Code	Internal diameter	Centre to centre
ST/94020	2013/2	50mm	85mm
ST/94022	2013/4	62mm	95mm
ST/94024	2013/5	69mm	95mm
ST/94026	2013/6	75mm	100mm
ST/94028	2013/8	82mm	103mm
ST/94030	2013/11	100mm	110mm
ST/94032	2013/14	125mm	125mm
ST/94034	2013/15	131mm	128mm
ST/94036	2013/17	150mm	135mm



ST/94020

Funnel, filter, holder, PP

Kartell

Single or double holder for funnels as small as 25mm (with adapter supplied) and up to 170mm. Without adapter can be used with powder funnels with stem diameter not exceeding 35mm. Two olives provided for use with rods from 8mm to 14mm.

Cat. No.	Code	Places	Funnel size, mm
FH/46222	268	1	25 - 170
FH/46224	269	2	25 - 170



FH/4622

Funnel, separating, holder, PP

Kartell

Horseshoe funnel holder especially for separating funnels. Supplied with two 'olives' so will fit supporting rod diameters 8 to 14mm. Cannot be used with funnels below 75mm diameter.

Cat. No.	Code	Description	Funnel size, ml
FH/46220	970	Funnel holder	125 - 1000



FH/46220

Clamp, burette

Kartell

A one-piece clamp which holds 50ml burettes firmly to a laboratory stand. Can be fitted and adjusted with one hand.

Cat. No.	Pk Qty
ST/44204	3



ST/44204