



SC/1620

Clamps, burette, polypropylene and rubber

- Spring loaded polypropylene clamps that are corrosion resistant
- Allow clear view of graduations along the full length of the burette.
- Rubber grips prevent slipping.
- A conical compression piece on the central boss provides a firm grip on rods 11 – 14mm in diameter

Cat. No.	Code	Description
SC/16202	CWH002	Single
SC/16204	CWH004	Double



ST/44202

Clamp, burette, double

High strength alloy casting with plastic jaws. Suitable for all burettes up to 100ml capacity.

Cat. No.	Description
ST/44202	Burette clamp, double




ST/44302

Stands, tripod, triangular, stainless steel

Cat. No.	Height x sides, mm
ST/94616	200 x 125

Stand, tripod, triangular, iron rod

Cat. No.	Height x sides, mm
ST/44302	200 x 125




ST/4801

Stands, triangles, pipeclay

Pipeclay stems on stout wire.

Cat. No.	Sides, mm	Pk Qty
ST/48012	50	10
ST/48014	63	10
ST/48016	75	10



ST/4410

Gauzes, plain

Supplied in packs of 10.

Cat. No.	Sides, mm
ST/44102	125 x 125
ST/44104	150 x 150