



SO/151

Overhead stirrers accessories

Stands



Extra strong retort stands with H-pattern base for stability and robust support rod. Both versions include bosshead.

Technical Specifications	SS10/1	SS10/2
Base (w x d x h)	400 x 350 x 25mm	550 x 480 x 25mm
Rod (dia. x l)	16 x 700mm	25 x 850mm
Net weight	7.6kg	20.2kg


Cat. No.	Code	Description
SO/15108	SS10/1	Stand, general purpose
SO/15110	SS10/2	Stand, heavy duty

Stirrer paddles

Made from high grade stainless steel, there is a choice of five paddle heads and two lengths of 8mm diameter paddle rod. The heads screw on and off the rods, so they can be mixed and matched depending upon stirring requirements. This gives maximum versatility and value for money. For complete paddle, order rod plus head(s)

Technical Specifications					
Paddle rods	SS10/5	SS10/6			
Diameter	8mm	8mm			
Length	350mm	550mm			
Paddle heads	SS10/10	SS10/11	SS10/12	SS10/13	SS10/14
Total width	60mm	94mm	80mm	50mm	60mm

Cat. No.	Code	Description
Paddle rods		
SO/15212	SS10/5	Paddle rod, 350mm
SO/15214	SS10/6	Paddle rod, 550mm
Paddle heads		
SO/15316	SS10/10	Small paddle head
SO/15318	SS10/11	Large paddle head
SO/15320	SS10/12	Large propeller head (2 blade)
SO/15322	SS10/13	Small propeller head (4 blade)
SO/15324	SS10/14	Anchor paddle head



(Paddle heads shown fixed to rod)



SO/4620

Stirring paddles, PP



A variety of blades made of chemically resistant polypropylene and fixed to a steel rod covered with Polypropylene. Rods are 8mm diameter and 350mm long.

Cat. No.	Code	Description	Blade w x h
SO/46202	428	Two fixed blades	39 x 15 mm
SO/46204	430	'U' shape blade	65 x 78 mm
SO/46206	431	Two collapsible blades	57 mm fully extended
SO/46208	432	Two collapsible blades	98.5 mm fully extended