

Stoppers, translucent silicone

Hardness 46/55 IRHD, to BS2775:87.

Cat. No.	Bottom diam. (mm)	Top diam. (mm)	Length (mm)	Pk Qty
SR/76302	3	4.5	12	50
SR/76304	4	5.8	14	50
SR/76306	5	7	16	50
SR/76308	6	8	16	50
SR/76310	7	9	16	50
SR/76312	8	10.5	20	20
SR/76314	9	11.5	20	20
SR/76316	10	12.5	20	20
SR/76318	11	14	24	10
SR/76320	13	16	24	10
SR/76322	15	18	24	10
SR/76324	17	20.3	26	10
SR/76326	19	22.5	28	5
SR/76328	21	24.5	28	5
SR/76330	23	26.5	28	5
SR/76332	25	28.5	28	2
SR/76334	27	31	32	2
SR/76336	29	33	32	2
SR/76338	31	36	35	2
SR/76340	33	38.5	38.5	2
SR/76342	35	40.5	38.5	1
SR/76344	37	42.5	38.5	1
SR/76346	39	44.5	38.5	1
SR/76348	41	47	42	1
SR/76350	43	49	42	1
SR/76352	45	51	42	1
SR/76354	49	55	42	1
SR/76356	53	59.5	45.5	1
SR/76358	57	64	49	1
SR/76360	62	69	49	1
SR/76362	67	74	49	1
SR/76364	71	78	49	1
SR/76370	86	96	54	1
SR/76374	99	109	54	1
SR/76376	105	116.5	63.5	1



SR/763

Larger sizes also available and with either one or two holes, please ask for details

Notes
