

Swab rinse kits for environmental sampling



Enables easy monitoring of bioburden levels on surfaces and equipment.

Sterile sampling templates (10cm x 10cm) allows accurate microbial loading of surfaces to be calculated. Alginate tipped swabs dissolve fully to ensure complete delivery of microbial sample into the SRK solution.

Cat. No.	Code	Description	Case Qty
SW/13372	902C	Swab rinse kit with standard alginate tip	25
SW/13374	903C	Swab rinse kit with alginate tip swab	25
SW/13376	904C	Swab rinse kit with large tip swab	25
SW/13378	905C	Square sampling kit	10

Swabs, amies, plastic stem

Swabs, Amies, Gamma irradiated, plastic stem, sterile, 12mm diameter x 140mm long, supplied in a polypropylene test tube, individually wrapped, peel pack of 100, expiry 560 days from date of manufacture..

Cat. No.	Code	Description
SW/90002	SV0021K	Swab Modified Amies Medium Transport
SW/90004	SV0023K	Swab Amies Medium Transport , with Charcoal

Swabs, disposable, dry tip

Sterilised by ethylene oxide. Individually wrapped.
For use when samples do not have to be transported.

Cat. No.	Code	Description	Case Qty
SW/23402	TS/8-A	Woodstick swab with a cotton tip	1000

Swabs, disposable, non sterile

Dimensions: 150 x 2.5 mm (tip \pm 5 mm).

Cat. No.	Code	Description	Case Qty
SW/23422	TS/10-NS	Woodstick swabs with cotton tip, 100 swabs sealed in a mini grip bag	5000

Swabs, disposable, sterile in shockproof tube

Sterilised by ethylene oxide. Individually wrapped in shockproof polypropylene round bottom tubes sealed by the label. Tube dimensions 13 x 165 mm (dia x length).

Cat. No.	Code	Description	Case Qty
SW/23532	TS/6-A	Woodstick swabs with cotton tip	500
SW/23534	TS/19-A	Polystyrene shafted swab with viscose tip	500

Swabs, disposable, sterile Viscose tip

Individually wrapped in a peel pouch, Gamma irradiated, flexible Polypropylene shaft with viscose tip providing a better sampling material than cotton.

Cat. No.	Description	Case Qty
SW/23406	Swab, polypropylene flexible shaft with a Viscose tip, sterile	1000

Wooden Applicators

Wooden applicators 150mm long packed in cases of 50 x 1000..


Cat. No.	Description
SW/43702	Wooden applicators



SW/133



SW/90002



SW/23532



SW/23534



SW/23406