




BW/65704



BW/65010



BW/6512



BW/6513



BW/85268

Baths, accessories

Lids

Recommended for use with temperatures above 60°C

Clear polycarbonate, gabled

Cat. No.	Code	Description
BW/65704	AQL5	To fit JBA5 bath
BW/65706	AQL12	To fit JBA12 bath
BW/65708	AQL18	To fit JBA18 bath

Stainless steel, gabled

Cat. No.	Code	Description
BW/65010	LU6	To fit JBA5 bath
BW/65012	LU14	To fit JBA12 bath
BW/65014	LU28	To fit JBA18 bath

For flat stainless steel lids with ring sets, please see BW/65026 to BW/65030 on page 42.

Test tube racks

For JB Aqua Plus and SUB Aqua Plus baths

Bath capacity	Number of J2 racks
5 litre	1
12 litre	2
18 litre	4
26 litre	4
Dual	1 + 2

Cat. No.	Code	Holes x Tube size
BW/65120	J2-10	84 x 10mmØ
BW/65122	J2-13	55 x 13mmØ
BW/65124	J2-16	36 x 16mmØ
BW/65126	J2-19	32 x 19mmØ
BW/65128	J2-25	18 x 25mmØ
BW/65130	J2-30	12 x 30mmØ

Microcentrifuge tube racks

Cat. No.	Code	Holes x Tube size
BW/65132	J2-SE	105 x 0.5ml
BW/65134	J2-LE	65 x 1.5ml

Spheres, polypropylene

Provide a floating lid on water baths to minimize evaporation and allows easy immersion of flasks etc

Number of packs required to provide a single layer

Bath model	Quantity of packs
SAP2, SAP2S, SAP5 and SAP12	1
SAP18 and SAP26	2
SAP34	3
SAPD	2
SBB AQUA 5 Plus and SBB AQUA 12 Plus	1
SBB AQUA 18 Plus and SBB AQUA 26 Plus	2
JBN5 and JBN12	1
JBN18 and JBN26	2
JBA5 and JB12	1
JBA18	2

Cat. No.	Description	Pk Qty
BW/85268	Spheres PP 20mm dia.	200