

Test tubes

Test tubes, soda glass



Samco soda glass test tubes are suitable for general mixing and simple laboratory work. They will withstand moderate warming (e.g. in a water bath). They should not be placed in naked flames.

Light wall, with rim

Cat. No.	Size, mm	Wall, mm	Approx, volume, ml	Pk Qty
TT/42002	50 x 6	0.45	0.2	250
TT/42004	50 x 10	0.6	2.8	250
TT/42006	75 x 10	0.6	4.2	100
TT/42008	75 x 12	0.6	6.2	100
TT/42010	100 x 12	0.6	8.3	100
TT/42012	100 x 16	0.6	15.5	100
TT/42014	125 x 16	0.6	19.4	100
TT/42016	125 x 18	0.6	25.0	100
TT/42018	150 x 16	0.6	23.3	100
TT/42020	150 x 18	0.6	30.0	100
TT/42022	150 x 24	0.8	53.2	50



TT/420

Medium wall, rimless, bacteriological quality

Cat. No.	Size, mm	Wall, mm	Approx, volume, ml	Pk Qty
TT/42104	50 x 10	1.0	2.3	250
TT/42106	75 x 10	1.0	3.4	100
TT/42108	75 x 12	1.0	5.4	100
TT/42110	100 x 12	1.0	7.1	100
TT/42112	125 x 16	1.0	17.4	100
TT/42118	150 x 16	1.0	20.8	100
TT/42120	150 x 18	1.0	27.2	100
TT/42122	150 x 24	1.2	49.5	50

Durham tubes



Soda glass, in packs of 250.

Cat. No.	Size, mm	Wall, mm	Approx, volume, ml
TT/42132	25 x 6.5	0.65	0.5
TT/42134	30 x 6.5	0.65	0.6
TT/42136	35 x 8	0.65	1.2
TT/42138	50 x 7.5	0.65	1.4

Test tubes, Samco '33' borosilicate glass



Samco borosilicate glass test tubes are suitable for all usual laboratory applications. They will withstand temperature changes up to and including boiling of samples; they are also highly chemically resistant.

Medium wall, with rim

Cat. No.	Size, mm	Wall, mm	Approx, volume, ml	Pk Qty
TT/42206	75 x 10	1.0	3.4	100
TT/42208	75 x 12	1.0	5.4	100
TT/42210	100 x 12	1.0	7.1	100
TT/42212	100 x 16	1.2	13.1	100
TT/42214	125 x 16	1.2	16.4	100
TT/42218	150 x 16	1.2	19.7	100
TT/42220	150 x 18	1.2	25.9	100
TT/42222	150 x 24	1.2	49.5	50

Test tubes, borosilicate glass, medium wall, rimless

Cat. No.	Size, mm	Wall, mm	Approx, volume, ml	Pk Qty
TT/42306	75 x 10	1.0	4.2	100
TT/42308	75 x 12	1.0	6.2	100
TT/42310	100 x 12	1.0	8.3	100
TT/42312	100 x 16	1.2	13.1	100
TT/42314	125 x 16	1.2	16.4	100
TT/42318	150 x 16	1.2	19.7	100
TT/42320	150 x 18	1.2	25.9	100
TT/42322	150 x 24	1.2	49.5	50



TT/423