



TT/124

Test tubes, Pyrex® borosilicate glass

PYREX®

High quality test tubes designed and constructed for optimum performance.

- Allows use at high temperatures
- Excellent resistance to chemical attack
- Heavy wall tubes have the additional advantage of lower replacement costs & improved user safety.

Medium wall with rim

Cat. No.	Size, mm	Wall, mm	Approx. volume, ml	Pk Qty
TT/12406	10 x 75	1.0	4.2	100
TT/12408	12 x 75	1.0	6.2	100
TT/12410	12 x 100	1.0	8.3	100
TT/12412	16 x 100	1.2	13.1	100
TT/12414	16 x 125	1.2	16.4	100
TT/12418	16 x 150	1.2	19.7	100
TT/12420	18 x 150	1.2	25.9	100
TT/12422	24 x 150	1.2	42.5	100

Medium wall, rimless


Cat. No.	Size, mm	Wall, mm	Approx. volume, ml	Pk Qty
TT/12506	10 x 75	1.0	4.2	100
TT/12508	12 x 75	1.0	6.2	100
TT/12510	12 x 100	1.0	8.3	100
TT/12512	16 x 100	1.2	13.1	100
TT/12514	16 x 125	1.2	16.4	100
TT/12518	16 x 150	1.2	19.7	100
TT/12520	18 x 150	1.2	25.9	100
TT/12522	24 x 150	1.2	42.5	100

Heavy wall with rim

Cat. No.	Size, mm	Wall, mm	Approx. volume, ml	Pk Qty
TT/12614	16 x 125	1.8	13.6	100
TT/12618	16 x 150	1.8	16.4	100
TT/12620	18 x 150	1.8	22.0	100
TT/12622	24 x 150	1.8	44.2	100

Heavy wall, rimless

Cat. No.	Size, mm	Wall, mm	Approx. volume, ml	Pk Qty
TT/12714	16 x 125	1.8	13.6	100
TT/12718	16 x 150	1.8	16.4	100
TT/12722	24 x 150	1.8	44.2	100



TT/120

Tubes, culture, Pyrex® borosilicate

PYREX®

Disposable rimless glass culture tubes. Inner pack 250. For full range please see page 462.

Cat. No.	Size, mm	Approx. volume, ml	Pk Qty
TC/12002	10 x 75	4.0	1000
TC/12004	12 x 75	6.0	1000
TC/12006	13 x 100	10.0	1000
TC/12008	16 x 100	11.0	1000
TC/12010	16 x 125	15.0	1000
TC/12012	16 x 150	19.0	1000
TC/12014	16 x 150	23.0	1000
TC/12016	18 x 150	28.5	500
TC/12018	20 x 150	36.0	500