

## Thermometers

Eco-friendly thermometers filled with various coloured non-toxic expansion liquid as a substitute for general purpose mercury filled thermometers. Fitted with a non-roll end cap with an integral ring to allow suspension of the thermometer in a sample. Supplied in a protective PVC casing.

Cat. No.	Code	Range, °C	Division, °C	Immersion	Length, mm	Pk Qty
<b>Green spirit</b>						
TH/33250	46/700/0	-35 to +50	1	76mm	305	5
TH/33252	46/700/1	-35 to +50	1	Total	305	5
TH/33254	46/701/0	-10 to +50	0.5	76mm	305	5
TH/33256	46/701/1	-10 to +50	0.5	Total	305	5
TH/33258	46/703/0	-10 to +110	1	76mm	305	5
TH/33260	46/703/1	-10 to +110	1	Total	305	5
TH/33262	46/706/0	-10 to +150	1	76mm	305	5
TH/33264	46/706/1	-10 to +150	1	Total	305	5
<b>Red spirit</b>						
TH/33230	44/815/0	-50 to 50	1	Total	300	Each
TH/33232	44/814/0	-100 to 50	1	Total	300	Each
TH/33202	44/430/0	-10 to 110	1	76mm	150	10
TH/33204	44/423/0	-10 to 110	1	76mm	205	5
TH/33206	44/401/0	-10 to 50	0.5	76mm	305	5
TH/33208	44/403/0	-10 to 110	1	76mm	305	5
TH/33210	44/406/0	0 to 150	1	76mm	305	5
<b>Blue spirit</b>						
TH/33312	44/830/7	-10 to 110	1	Total	155	Each
TH/33314	44/808/8	-10 to 50	0.5	Total	305	Each
TH/33316	44/809/8	-35 to 50	1	Total	305	Each
TH/33318	44/810/8	-20 to 110	1	Total	305	Each
TH/33320	44/811/8	-20 to 150	1	Total	305	Each
TH/33322	44/812/8	-10 to 210	1	Total	405	Each
TH/33342	44/800/8	-10 to 50	0.5	76mm	305	Each
TH/33344	44/801/8	-35 to 50	1	76mm	305	Each
TH/33346	44/802/8	-20 to 110	1	76mm	305	Each
TH/33348	44/803/8	-20 to 150	1	76mm	305	Each
TH/33350	44/804/8	-10 to 260	1	76mm	405	Each
<b>Brown spirit</b>						
TH/33352	44/815/0	-10 to 300	1	Total	300	Each
TH/33354	44/814/0	-10 to 330	2	Total	300	Each

\*\*Other ranges available, please enquire



Green spirit



Red spirit

## Thermometer case, slot front, brass

For thermometers 5-7mm in diameter.

Cat. No.	Length accommodated, mm
TH/33402	155
TH/33404	305
TH/33406	355
TH/33408	405



TH/3340