



TH/33014

## Thermometer, fridge freezer

White horizontal fridge/freezer thermometer with blue temperature scale indicating recommended zones.

Cat. No.	Description
TH/33014	Fridge/freezer thermometer



TH/95150

## Thermometer, fridge freezer with audible alarm

Digital fridge freezer alarm with a water-resistant probe with a 1-meter cable. Records maximum and minimum temperatures from -24.9°C to 69.9°C.

- Programmable high/low audible alarm
- Twin sensors measure external and internal temperature
- Audible temperature alarm
- Accuracy of  $\pm 1^\circ\text{C}$
- Large easy to read LCD

Cat. No.	Code	Description
TH/95150	810-125	Thermometer, fridge freezer, -24.9 to 69.9°C



TH/95160



TH/95162


## Thermolite2 digital thermometer

Housed in a high impact ABS plastic case with a built-in stainless-steel probe. The unit features a CalCheck  $0.0^\circ\text{C}$  ( $\pm 0.1^\circ\text{C}$ ) function that allows the user to verify the accuracy of the thermometer. The addition of a silicone boot, available separately, increases the water resistance to IP65 rating. Supplied with a free traceable certificate of calibration.

- Temperature range  $-39.9$  to  $149.9^\circ\text{C}$ , resolution  $0.1^\circ\text{C}$ , accuracy of  $\pm 0.5^\circ\text{C}$  ( $-39.9$  to  $124.9^\circ\text{C}$ )
- Low battery indication
- High impact ABS plastic case
- Available in six colours, suitable for different applications
- Tip diameter  $3.3\text{mm} \times 80\text{mm}$  long
- Clear easy to read display

Cat. No.	Code	Description
TH/95160	226-112	Thermolite2 pocket digital thermometer, White* with built-in stainless-steel probe
TH/95162	830-280	Boot, silicone for Thermolite2 for water resistance to IP65 rating

\*Also available in Grey, Yellow, Green, Red or Blue. Please specify colour required, White is automatically supplied if no colour is specified.



TH/35012

## Thermometer, Key pocket



A K-type pocket thermometer, which is available with a selection of interchangeable probes. Robust and easy to operate, the Key is supplied complete with HI-98517-13 general purpose/penetration probe and batteries.

Technical Specifications	
Range	-40 to $550^\circ\text{C}$
Resolution	$1^\circ\text{C}$
Accuracy	$\pm 2^\circ\text{C}$
Dims/weight	$175 \times 41 \times 23\text{mm}/78\text{g}$

Cat. No.	Code	Description
TH/35012	HI-98517	Keypocket thermometer
<b>Accessories</b>		
TH/35018	HI-98517-13	Penetration probe 130mm
TH/35020	HI-98517-15	General purpose probe 170mm
TH/35022	HI-98517-30	General purpose probe 320mm
TH/35024	CERT-CAL	UKAS traceable calibration certificate