



TW/95202

Stopwatches, digital

Housed in an impact-resistant ABS case these stopwatches are easy to use for timing single events.

ATM-136 7.2mm Digit Stopwatch

- 1/100 sec. precision with lap counter (for first 30 minutes) and 1 second thereafter
- Calendar & time (12/24-hour format)
- Daily alarm
- Hourly chime
- Big digit display (digit height 7.2mm)
- Case size: 64 x 53 x 19mm
- Supplied with battery and neck cord

ATM-308 12mm Digit Stopwatch

- 1/100 sec. precision with lap counter (for first 30 minutes) and 1 second thereafter
- Calendar & time (12/24-hr format)
- Daily alarm
- Hourly chime
- Water and shock resistant structure
- Case size: 72 x 64 x 20mm
- Supplied with battery and neck cord

Cat. No.	Code	Description
TW/95202	ATM-136	7.2mm Digit Stopwatch
TW/95204	ATM-308	12mm Digit Stopwatch
TW/95208	ACAL-TM	Calibration certificate for stopwatch



TW/65014

Timer, count down

Timer counts down from 60 minutes with alarm, 60 x 60 x 35mm. Accepts R6 battery (supplied).

Cat. No.	Description
TW/65014	Timer, count down

Timer, count up and down, digital

A simple to use count-up or count-down timer. It is ideal for laboratories, education and industry.

- 99min's 99sec's x 1 second
- Large 20mm LCD display
- Count-down mode with Alarm on expiry
- Count-down can be stopped and restarted at any point
- Memory to store last count-down time setting
- Repeat function for Count-down times
- Cumulative count-up function – count up can be paused and re-started without losing the cumulative time count
- Impact-resistant ABS case with flip-out desk stand and rear magnet



TW/95316

Cat. No.	Description
TW/95316	Timer, count up/down

Timer, count up and down, digital, three channel

Counts up and down to/from 24 hours with alarm at end of count down period, three timing channels, clock with stopwatch, memory, 10mm digits, stand and magnet, 93x70mm, batteries 2 x LR44.



TW/65122

Cat. No.	Description
TW/65122	Timer, count up/down