



TC/160



TC/12020



TC/12032



TC/12046

## Tissue Culture



### Culture tubes, disposable plastic

Suitable for suspension culture procedures

- Available in two sizes: 12x75mm (5ml) and 17 x 100mm (13ml)
- Feature a two position cap, leave loose for gaseous exchange or push down tight to seal
- Available in polypropylene or polystyrene with a round base
- Individually or bulk packed
- All tubes feature a white marker spot for easy labelling and printed graduations
- Tubes are gamma irradiated for ensured sterility

Cat. No.	Product Code	Description	Tube size, mm	Approx volume, ml	Inner Pack	Case Qty
<b>Polystyrene tubes</b>						
TC/16002	75CTS1	5ml culture tube	12x75	4	1x100	500
TC/16004	75CTS	5ml culture tube	12x75	4	25	500
TC/16006	100CTS1	13ml culture tube	17x100	12	1x100	500
TC/16008	100CTS	13ml culture tube	17x100	12	25	500
<b>Polypropylene tubes</b>						
TC/16010	75CTP1	5ml culture tube	12x75	4	1x100	500
TC/16012	75CTP	5ml culture tube	12x75	4	25	500
TC/16014	100CTP1	13ml culture tube	17x100	12	1x100	500
TC/16016	100CTP	13ml culture tube	17x100	12	25	500

### Culture Tubes, PYREX® - disposable glass

**PYREX®**

Manufactured from chemically inert Pyrex® glass. Tubes and caps can withstand autoclaving at 121°C.

Cat. No.	Product Code	Dimensions, mm	Approx volume, ml	Case Qty
<b>Glass Culture Tube, Rimless, Round Bottom</b>				
TC/12002	99445-10	10x75	4	1000
TC/12004	99445-12	12x75	6	1000
TC/12006	99445-13	13x100	10	1000
TC/12008	99445-15	15x85	11	1000
TC/12010	99445-16	16x100	15	1000
TC/12012	99445-16X	16x125	19	1000
TC/12014	99445-16XX	16x150	23	1000
TC/12016	99445-18	18x150	28.5	500
TC/12018	99445-20	20x150	36	500
<b>Glass Culture Tube, Screw Neck, Round Bottom</b>				
TC/12020	99449-13	13x100	8.5	1000
TC/12022	99449-16	16x100	15	1000
TC/12024	99449016X	16x125	19	1000
TC/12026	99449-16XX	16x150	20	1000
TC/12028	99449-20	20x125	30	500
TC/12030	99449-20X	20x150	34	500
<b>Glass Culture Tube, Screw Neck, Flat Bottom</b>				
TC/12032	99448-16	16x125	17	1000
TC/12034	99448-19	19.5x145	29.5	500
<b>Screw cap with rubber liner</b>				
TC/12040	99999-13	13		1000
TC/12042	99999-15	16		1000
TC/12044	99999-18	20		1000
<b>Screw cap with PTFE faced liner</b>				
TC/12046	9998-13	13		288
TC/12048	9998-15	16		288
TC/12050	9998-18	20		192