



TC/463 Series

### Stripette, individual, paper/plastic wrap, Corning

Paper/plastic wrap allows easy opening by either pop-through or peel-apart techniques. It is also ideal for sterile tissue culture applications, especially while wearing latex gloves, as the wrap reduces static cling. Colour coding on paper allows easy size identification.

Cat. No.	Code	Capacity, ml	Grad., ml	Negative Grad., ml	Colour	Inner Packs	Case Qty
TC/46310	4485	1	1/100	0.2	Yellow	50	1000
TC/46312	4486	2	1/100	0.2	Green	50	1000
TC/46318	4487	5	1/10	2.5	Blue	50	200
TC/46328	4488	10	1/10	3	Orange	50	200
TC/46334	4489	25	2/10	10.0	Red	25	200
TC/46340	4490	50	1/1	10.0	Purple	25	100
TC/46346	4491	100	1/1	n/a	Aqua	10	100



TC/466 Series

### Stripette, bulk pack, Corning

Bulk pack pipettes are ideal for large scale sterile and non-sterile liquid handling applications.

Cat. No.	Code	Capacity, ml	Grad., ml	Negative Grad. ml	Colour	Inner Packs	Case Qty
TC/46610	4010	1	1/100	0.2	Yellow	50/bag	1000
TC/46612	4020	2	1/100	0.2	Green	25/bag	1000
TC/46618	4050	5	1/10	2.5	Blue	25/bag	500
TC/46628	4100	10	1/10	3.0	Orange	25/bag	500
TC/46634	4250	25	2/10	10.0	Red	20/bag	200
TC/46640	4500	50	1/2	10.0	Purple	20/bag	100




TC/4680

### Polypropylene container, Corning

Corning flexible polypropylene container with snap-on polythene lid serves as a beaker or storage container. Graduated in both millilitres and ounces. Certified non-pyrogenic. Supplied sterile.

Cat. No.	Code	Description	Capacity, ml	Inner Packs	Case Qty
TC/46802	430179	Container & Lid	250ml	1	100
TC/46804	430180	Container only	250ml	20	500
TC/46806	430181	Lid only	n/a	20	500



TC/46846

### Graduated cylinder, polystyrene, Corning

Optically clear polystyrene graduated cylinders are designed for sterile, accurate dispensing of culture media or other biological fluids. A polyethylene dust cover is included. Supplied sterile.

Cat. No.	Code	Capacity, ml	Grads., ml	Inner Packs	Case Qty
TC/46846	430182	100	1	1	50



TC/468

### Erlenmeyer flask, polycarbonate Corning

Made from optically clear polycarbonate. Ideal for shaker culture applications. Sterilized by gamma radiation and certified non pyrogenic.

Cat. No.	Code	Capacity, ml	Grads., ml	Neck diam.	Inner Packs	Case Qty
<b>Plug seal cap</b>						
TC/46848	430421	125	25	26 mm	1	50
TC/46856	430183	250	25	31 mm	1	50
TC/46866	430422	500	50	43 mm	1	25
TC/46876	431146	1000	50	43 mm	1	25
<b>Vented cap</b>						
TC/46948	431143	125	25	26 mm	1	50
TC/46956	431144	250	25	31 mm	1	50
TC/46966	431145	500	50	43 mm	1	25
TC/46976	431147	1000	50	43 mm	1	25