

Trays, tidy, PVC

Kartell

Very useful for small objects - rulers, erasers, pencils etc. in desk drawers. Central compartment of five compartment tray can be used for pipettes. Sturdy.

Cat. No.	Code	Compartments	Dimensions (WxDxH), mm
TY/46352	953	12 (90 x 90mm)	403 x 303 x 63
TY/46354	954	4 x (100 x 185mm) 1 x (100 x 385mm)	404 x 304 x 64



TY/4635

Utility carrier

NALGENE

Rigid utility carrier with sturdy handle and moulded in centre division. Useful for carrying a multitude of lab items. Available in HDPE and Polypropylene. The Polypropylene version can be autoclaved and meets OSHA Standard 29 CFR Part 1910.1030 for use as protection against blood borne pathogens. Overall dimensions (including handle), L x W x H, 378 x 241 x 114mm. Compartment dimensions, L x W x H, 340 x 89 x 111mm.

Cat. No.	Code	Description
BC/36102	7130-0010	Utility carrier HDPE
BC/36204	7134-0010	Utility carrier PP



BC/3610

Trolleys

For gas cylinder trolleys please see page 220.

Trolley, shelf

General purpose trolley with a grey mild steel frame and removable aluminium shelves. All swivel castors. Overall dimensions; 800 x 550 x 915mm high. Height to top working surface; 770mm. Tray size: 735 x 530mm with 40mm upstand

Cat. No.	Description
TR/94080	Trolley, 2 Tier, clearance between shelves: 560mm
TR/94082	Trolley, 3 Tier, clearance between shelves: 260mm



TR/94080

Trolley, multi-purpose

Can be used to transport a wide range of items within schools, hospitals or laboratories. Powder coated mild steel frame with 2 removable plastic coated wire baskets. All swivel castors. Basket size: 500 x 388 x 170mm (LxWxD). Overall dimensions: 625 x 442 x 1020mm (LxWxH)

Cat. No.	Description
TR/94090	Multi-purpose trolley with 2 baskets



TR/94090