

Tubing

Tubing, soda glass

Light wall, lengths of 500mm.

Cat. No.	O.D., mm	Pk Qty
TU/02002	4	30
TU/02004	5	30
TU/02006	6	30
TU/02008	7	30
TU/02010	8	30
TU/02012	10	30




TU/020

Tubing, borosilicate glass

Medium wall (1.0mm), lengths of 500mm.

Cat. No.	O.D., mm	Pk Qty
TU/02104	6	30
TU/02106	7	30
TU/02108	8	30
TU/02110	10	30



TU/021 Series

Tubing, borosilicate glass, capillary

Lengths of 500mm.

Cat. No.	O.D., mm	Bore, mm	Pk Qty
TU/02204	5	0.5	10
TU/02206	6	0.8	10
TU/02210	7	1.2	10



TU/4652

Tubing, PVC, clear, non toxic, flexible, reinforced

Cat. No.	Bore, mm	Metres	Wall thickness max
TU/46524	6.5	30	3
TU/46526	8	30	3
TU/46528	10	30	3
TU/46530	12.5	30	3
TU/46532	16	30	3
TU/46534	19.0	30	4



TU/765

Flame retardant bunsen tubing, black nitrile with reinforced ends

Supplied Individually.

Cat. No.	Code	Length, mm
TU/76502	CBT/450	450
TU/76504	CBT/600	600
TU/76506	CBT/750	750
TU/76508	CBT/900	900