




TU/761

Tubing, red rubber, B.S. 2775:1987

Cat. No.	Code	Bore, mm	Wall,mm	Metres
TU/76102	RT/N03	3.0	0.75	10
TU/76104	RT/N05	5.0	1.5	10
TU/76106	RT/N06.5	6.5	1.5	10
TU/76108	RT/N08	8.0	2	10
TU/76110	RT/N10	10.0	2	10
TU/76112	RT/N12.5	12.5	2.25	10
TU/76114	RT/N16	16.0	3.25	10
TU/76116	RT/N20	20.0	3.25	10
TU/76118	RT/N25	25.0	3.25	10



TU/762

Tubing, red rubber, heavy wall, pressure, BS2775:1987

Cat. No.	Code	Bore, mm	Wall,mm	Metres
TU/76220	RT/H03	3	3	10
TU/76222	RT/H05	5	4	10
TU/76224	RT/H06.5	6.5	5	10
TU/76226	RT/H08	8	6	10
TU/76228	RT/H10	10	7	10
TU/76230	RT/H12.5	12.5	8.5	10



TU/763

Tubing, Neoprene, black BS2775

Cat. No.	Code	Bore, mm	Wall,mm	Metres
TU/76342	NT/N03	3.0	0.75	10
TU/76344	NT/N05	5.0	1.5	10
TU/76346	NT/N06.5	6.5	1.5	10
TU/76348	NT/N08	8.0	2.0	10
TU/76350	NT/N10	10.0	2.0	10
TU/76352	NT/N12.5	12.5	2.0	10

Tubing, translucent silicone

Platinum cured, 60 S.H., medical grade, non-toxic, autoclavable by steam at 130°C. Coil of 10m.

Cat. No.	Bore, mm	Wall, mm	O.D. (mm)
TU/76462	3	0.75	4.5
TU/76464	5	1.5	8
TU/76466	6.5	1.5	9.5
TU/76467	8	2	12
TU/76468	10	2	14
TU/76470	12.5	2.25	17
TU/76472	16	3.25	22.5