

## Tubing, Peristaltic pump




Micro bore tubing for peristaltic pumps has a reputation for reliable performance. Tubing has been designed to fit many peristaltic pumps and instrumentation which include; Pumps: Ismatec™, Gilson®, Watson-Marlow® and Skalar to name a few.

Instruments: Perkin Elmer™, Thermo® ICP instrumentation, flame photometers and many others. Ideal for research and routine use. A range of raw materials are listed which offer varying degrees of chemical resistance.

### Accurated PVC tubing

Flow rated and heat treated for increased life. They are suitable for general laboratory use and are resistant to most routine chemicals. Transparency allows for convenient flow monitoring and PVC tubing bends easily, grips connectors firmly and simplifies pump set up. Pack of 12.

| Cat. No. | Bridge colour   | Nominal OD, mm | Flow cc/min<br>(typical average) |
|----------|-----------------|----------------|----------------------------------|
| TU/51002 | Orange / Red    | 2.019          | 0.03                             |
| TU/51004 | Orange / Blue   | 2.083          | 0.05                             |
| TU/51006 | Orange / Green  | 2.210          | 0.10                             |
| TU/51008 | Orange / Yellow | 2.337          | 0.16                             |
| TU/51010 | Orange / White  | 2.464          | 0.23                             |
| TU/51012 | Black / Black   | 2.438          | 0.32                             |
| TU/51014 | Orange / Orange | 2.565          | 0.42                             |
| TU/51016 | White / White   | 2.692          | 0.60                             |
| TU/51018 | Red / Red       | 2.819          | 0.80                             |
| TU/51020 | Grey / Grey     | 2.972          | 1.00                             |
| TU/51022 | Yellow / Yellow | 3.099          | 1.20                             |
| TU/51024 | Blue / Blue     | 3.327          | 1.60                             |
| TU/51026 | Green / Green   | 3.531          | 2.00                             |
| TU/51028 | Purple / Purple | 3.734          | 2.50                             |
| TU/51030 | Purple / Black  | 3.962          | 2.90                             |
| TU/51032 | Purple / Orange | 4.216          | 3.40                             |
| TU/51034 | Purple / White  | 4.470          | 3.90                             |
| TU/51036 | Blue / Yellow   | 3.200          | 1.40                             |



TU/510

## Tubing, P.V.C. clear, very flexible

Can be used with peristaltic pumps.

| Cat. No. | Bore, mm | Wall, mm | Metres |
|----------|----------|----------|--------|
| TU/46502 | 5        | 1.5      | 30     |
| TU/46504 | 6.5      | 1.5      | 30     |
| TU/46505 | 8        | 2        | 30     |
| TU/46506 | 10       | 2        | 30     |
| TU/46508 | 12.5     | 2.25     | 10     |
| TU/46510 | 16       | 3.25     | 10     |
| TU/46512 | 20       | 3.25     | 10     |



TU/4650

## Tubing, pump, standard PVC



Supplied in 15m roll.

| Cat. No. | Bore, mm | Nominal OD, mm |
|----------|----------|----------------|
| TU/51102 | 0.254    | 2.083          |
| TU/51104 | 0.381    | 2.210          |
| TU/51106 | 0.508    | 2.337          |
| TU/51108 | 0.635    | 2.464          |
| TU/51110 | 0.762    | 2.438          |
| TU/51112 | 0.889    | 2.565          |
| TU/51114 | 1.016    | 2.692          |
| TU/51116 | 1.143    | 2.819          |
| TU/51118 | 1.295    | 2.972          |
| TU/51120 | 1.422    | 3.099          |
| TU/51122 | 1.651    | 3.327          |
| TU/51124 | 1.854    | 3.531          |
| TU/51126 | 2.057    | 3.734          |
| TU/51128 | 2.286    | 3.962          |
| TU/51130 | 2.540    | 4.216          |
| TU/51132 | 2.794    | 4.470          |
| TU/51134 | 1.524    | 3.200          |
| TU/51136 | 3.175    | 4.851          |



TU/511