

## Tubing connectors, 'T' pattern

Polypropylene, with rigid tapered arms. Autoclavable. Good chemical resistance.

*Azlon*

Cat. No.	Code	O.D. mm	I.D. mm	Pk Qty
TH/16102	CWH204	4 - 7	2	10
TH/16104	CWH206	6 - 7	4	10
TH/16106	CWH208	8 - 9	7	10
TH/16108	CWH210	10 - 12	9	10



TH/1610

## Tubing connectors, 'Y' pattern

Polypropylene, with rigid tapered arms. Autoclavable. Good chemical resistance.

*Azlon*

Cat. No.	Code	O.D. mm	Pk Qty
TH/16202	CWH304	4 - 7	10
TH/16204	CWH306	6 - 7	10
TH/16206	CWH308	8 - 9	10
TH/16208	CWH310	10 - 12	10



TH/1620

## Tubing connectors, straight pattern

Polypropylene with stepped tapered arms. Autoclavable. Good chemical resistance.

*Azlon*

Cat. No.	Code	O.D. mm	Pk Qty
TH/16302	CWH600	4-7 and 4-7	10
TH/16304	CWH602	4-7 and 7-12	10
TH/16306	CWH606	4-7 and 10-15	10
TH/16308	CWH604	7-12 and 7-12	10
TH/16310	CWH608	7-12 and 10-15	10
TH/16312	CWH610	10-15 and 10-15	10



TH/163