

CPX Series Ultrasonic Baths (continued)



UB/952



UB/953



UB/954



UB/955



Accessories

Technical Specifications

Model	1800	2800	3800	5800	8800
Tank capacity, litres	1.9	2.8	5.7	9.5	20.8
Tank size (LxWxD), mm	150x140x100	230x140x100	290x150x150	290x240x150	495x290x150
Overall size, (LxWxD), mm	251x304x302	337x304x302	397x317x377	397x401x378	596x466x391
Weight, Kg	7	9	12	14	26
Frequency	40Hz	40Hz	40Hz	40Hz	40Hz
Drain tap	No	No	Yes	Yes	Yes

Cat. No. Code Description

Bransonic 1800, 1.9L capacity

UB/95122	CPX-952-136R	Bath ultrasonic Branson Model M1800-E, mechanical timer only
UB/95124	CPX-952-137R	Bath ultrasonic Branson Model M1800H-E, mechanical timer - heater on/ off, 55°C
UB/95128	CPX-952-139R	Bath ultrasonic Branson Model CPX1800-E, digital timer only
UB/95126	CPX-952-138R	Bath ultrasonic Branson Model CPX1800H-E, digital timer - adjustable heater 20 - 69°C
UB/95110	100-410-170	Tray, solid, stainless steel 110 x 110 x 70mm
UB/95112	000-019-129	Wire basket, stainless steel, 125 x 107 x 74mm
UB/95114	CPN-246-010	Beaker holder, plastic, 1 x 600ml

Bransonic 2800, 2.8L capacity

UB/95222	CPX-952-236R	Bath ultrasonic Branson Model M2800-E, mechanical timer only
UB/95224	CPX-952-237R	Bath ultrasonic Branson Model M2800H-E, mechanical timer - heater on/ off, 55°C
UB/95228	CPX-952-239R	Bath ultrasonic Branson Model CPX2800-E, digital timer only
UB/95226	CPX-952-238R	Bath ultrasonic Branson Model CPX2800H-E, digital timer - adjustable heater 20 - 69°C
UB/95210	100-410-172	Tray, solid, stainless steel 190 x 90 x 75mm
UB/95212	000-019-229	Wire basket, stainless steel, 125 x 107 x 74mm
UB/95214	CPN-246-011	Beaker holder, plastic, 2 x 250ml
UB/95216	CPN-246-015	Beaker holder, plastic, 2 x 600ml

Bransonic 3800, 5.7L capacity

UB/95322	CPX-952-336R	Bath ultrasonic Branson Model M3800-E, mechanical timer only
UB/95324	CPX-952-337R	Bath ultrasonic Branson Model M3800H-E, mechanical timer - heater on/ off, 55°C
UB/95328	CPX-952-339R	Bath ultrasonic Branson Model CPX3800-E, digital timer only
UB/95326	CPX-952-338R	Bath ultrasonic Branson Model CPX3800H-E, digital timer - adjustable heater 20 - 69°C
UB/95310	100-410-174	Tray, solid, stainless steel 220 x 110 x 100mm
UB/95312	000-019-329	Wire basket, stainless steel, 261 x 119 x 136mm
UB/95314	CPN-246-012	Beaker holder, plastic, 3 x 250ml
UB/95316	CPN-246-016	Beaker holder, plastic, 2 x 600ml

Bransonic 5800, 9.5L capacity

UB/95422	CPX-952-536R	Bath ultrasonic Branson Model M5800-E, mechanical timer only
UB/95424	CPX-952-537R	Bath ultrasonic Branson Model M5800H-E, mechanical timer - heater on/ off, 55°C
UB/95428	CPX-952-539R	Bath ultrasonic Branson Model CPX5800-E, digital timer only
UB/95426	CPX-952-538R	Bath ultrasonic Branson Model CPX5800H-E, digital timer - adjustable heater 20 - 69°C
UB/95410	100-410-176	Tray, solid, stainless steel 250 x 190 x 150mm
UB/95412	000-019-529	Wire basket, stainless steel, 269 x 203 x 138mm
UB/95414	CPN-246-013	Beaker holder, plastic, 4 x 600ml

Bransonic 8800, 20.8L capacity

UB/95522	CPX-952-836R	Bath ultrasonic Branson Model M8800-E, mechanical timer only
UB/95524	CPX-952-837R	Bath ultrasonic Branson Model M8800H-E, mechanical timer - heater on/ off, 55°C
UB/95528	CPX-952-839R	Bath ultrasonic Branson Model CPX8800-E, digital timer only
UB/95526	CPX-952-838R	Bath ultrasonic Branson Model CPX8800H-E, digital timer - adjustable heater 20 - 69°C
UB/95510	100-410-178	Tray, solid, stainless steel 450 x 240 x 150mm
UB/95512	000-019-829	Wire basket, stainless steel, 471 x 271 x 130mm
UB/95514	CPN-246-014	Beaker holder, stainless steel, 6 x 600ml

Accessories for all models

UB/95516	000-019-226	Glass beaker, 250ml
UB/95518	000-019-125	Glass beaker, 600ml
UB/95520	000-583-700	O ring for 250ml beaker
UB/95530	000-583-710	O ring for 600ml beaker

Either basket, tray or breaker holder should be selected for correct operation.