

## Vials

### Vials, screw neck, squat form



Neutral glass vials, squat form, without closure. Polypropylene caps must be ordered separately.

| Cat. No. | Nominal capacity, ml | Ht x dia mm | Cap reference | Cap code | Pk Qty |
|----------|----------------------|-------------|---------------|----------|--------|
| VI/42006 | 7.0                  | 42 x 20     | VC/42208      | D        | 264    |
| VI/42008 | 14.0                 | 58 x 23.0   | VC/42212      | F        | 190    |
| VI/42010 | 28.25                | 72 x 27.5   | VC/42214      | G        | 140    |




VI/421

### Vials, screw neck, tall form



Neutral glass vials, tall form, without closure. Polypropylene caps must be ordered separately.

| Cat. No. | Nominal capacity, ml | Ht x dia mm | Cap reference | Cap code | Pk Qty |
|----------|----------------------|-------------|---------------|----------|--------|
| VI/42102 | 1.75                 | 36 x 11     | VC/42202      | A        | 882    |
| VI/42104 | 3.5                  | 46 x 12.5   | VC/42204      | B        | 666    |
| VI/42106 | 7.0                  | 58 x 17     | VC/42206      | C        | 399    |
| VI/42108 | 10.5                 | 67 x 17     | VC/42206      | C        | 399    |
| VI/42110 | 14.0                 | 72 x 20     | VC/42208      | D        | 264    |
| VI/42112 | 21.25                | 86 x 21     | VC/42210      | E        | 231    |
| VI/42114 | 28.25                | 96 x 23     | VC/42212      | F        | 190    |



VC/422

### Vials, rolled rim



Neutral glass vials without closure. These vials will accept either a polythene push-in plug or a clip on closure which must be ordered separately.

| Cat. No.              | Nominal capacity, ml   | Ht x dia mm | Cap reference         | Cap code | Pk Qty |
|-----------------------|------------------------|-------------|-----------------------|----------|--------|
| VI/42230              | 7.0                    | 34 x 23     | VC/46102 & 12         | H        | 190    |
| VI/42232              | 10.5                   | 46 x 23     | VC/46102 & 12         | H        | 190    |
| VI/42234              | 14.0                   | 48 x 25.5   | VC/46104 & 14         | J        | 162    |
| VI/42236              | 21.25                  | 66 x 25.5   | VC/46104 & 14         | J        | 162    |
| VI/42238              | 28.25                  | 72 x 27.5   | VC/46106 & 16         | K        | 140    |
| <b>Caps For vials</b> |                        |             | <b>Vial reference</b> |          |        |
| VC/46102              | Clip on cap, polythene |             | VI/42230 & 232        | H        | 190    |
| VC/46104              | Clip on cap, polythene |             | VI/42234 & 236        | J        | 162    |
| VC/46106              | Clip on cap, polythene |             | VI/42238              | K        | 140    |
| VC/46112              | Push in cap, polythene |             | VI/42230 & 232        | H        | 190    |
| VC/46114              | Push in cap, polythene |             | VI/42234 & 236        | J        | 162    |
| VC/46116              | Push in cap, polythene |             | VI/42238              | K        | 140    |



VI/4223 and VC/461