




WG/9200

## Watch glasses

### Watch glasses, glass

Manufactured from soda glass with ground edges.

Cat. No.	Diameter mm	Pk Qty
WG/92002	40	10
WG/92004	50	10
WG/92006	60	10
WG/92008	70	10
WG/92010	80	10
WG/92012	90	10
WG/92014	100	10
WG/92016	130	10
WG/92018	150	10



WG/120

### Watch glasses, Pyrex

**PYREX®**

For evaporation work and covering heated beakers. Fire-polished edges to reduce chipping and breakage.

Cat. No.	Item Code	Nominal diameter mm
WG/12002	1060/04D	40
WG/12004	1060/06D	50
WG/12006	1060/09D	60
WG/12008	1060/10D	80
WG/12010	1060/14D	100
WG/12012	1060/18D	125
WG/12014	1060/20D	150
WG/12016	1060/21D	200
WG/12018	1060/24D	250



WG/462

### Watch glasses, PP

**Kartell**

Polypropylene watch glasses with moulded base for stability. Can withstand temperatures up to 120°C.

Cat. No.	Code	Diameter mm	Height mm
WG/46206	452	60	10
WG/46210	454	80	12
WG/46214	455	100	14
WG/46215	456	120	18

## Notes

---



---



---



---



---



---



---