



WP/44100 & WP/44152



WP/44154



WP/44200



WP/44204

PURELAB®, chorus Reservoirs



Select a storage reservoir designed to maintain the optimum purity of stored purified water in a choice of 15, 30, 60 and 100 liter capacities.

Cat. No.	Code	Description
WP/44150	LA757	Reservoir, 15 litres
WP/44152	LA758	Reservoir, 30 litres
WP/44154	LA759	Reservoir, 60 litres
WP/44156	LA760	Reservoir, 100 litres

Each reservoir supplied with LA773 Installation Kit & LC216 Composite Vent Filter

Consumables

Cat. No.	Code	Description
WP/44170	LC216	Composite vent filter
WP/44172	LC275	Purification pack
WP/44174	LC210	Lamp, UV

**It is recommended that the filter is changed every 6 months, the purification pack every 6 months and the UV lamp every 12 to 18 months

Accessories

Cat. No.	Code	Description
WP/44180	LA770	Wall Mounting Kit 15/30 L Reservoir
WP/44182	LA771	Wall Mounting Kit 60/100 L Reservoir

PURELAB®, CHORUS HALO dispensers



The modular nature of PURELAB Chorus means that your dispense solutions sit independently from the unit. You can even have the Halo Dispenser installed in an adjacent laboratory.

HALO	Dispense	Advanced	Flexible
Fixed handset	✓	✓	
Flexible handset			✓
Variable flow rate dispense	✓	✓	✓
Drop-by-drop control	✓	✓	✓
Locked dispense	✓	✓	✓
Purity monitoring to point-of-use		✓	✓
Auto volume dispense		✓	✓

Cat. No.	Code	Description
WP/44200	LA754	Dispenser, Halo
WP/44202	LA755	Dispenser, Advanced, Halo
WP/44204	LA756	Dispenser, Flexible, Halo

Each dispenser is supplied with GUID40007 or GUID40006 Quick Reference Guide, INST40011 or INST40013 Quick Start Guide, LA774 Installation Kit (plus for LA756: MANU40002 Operator manual)

Accessories

Cat. No.	Code	Description
WP/44210	LA768	Mounting Kit, remote, for all Models

Consumables

Cat. No.	Code	Description
WP/44220	LC145	Filter, point of use
WP/44222	LC197	Biofilter point of use

It is recommended that the filters are changed every 3 months