



BE/162

Beakers, plastic

Beakers, tapered style, moulded graduation, polypropylene

- Tapered stacking beakers in high clarity PP.
- Designed to comply with ISO7056.
- High clarity PP is autoclavable at 121°C.
- Durable and economical with distinct moulded graduations.
- Precise non-drip pouring spouts.
- S.P.I. recycling code and date of manufacture imprints.

Cat. No.	Code	Capacity, ml	Sub divs, ml	Pk Qty
BE/16234	BPM0025P	25	5	10
BE/16240	BPM0050P	50	10	10
BE/16246	BPM0100P	100	20	10
BE/16256	BPM0250P	250	50	10
BE/16262	BPM0400P	400	50	10
BE/16266	BPM0500P	500	50	10
BE/16270	BPM0600P	600	50	5
BE/16276	BPM1000P	1000	100	5
BE/16280	BPM2000P	2000	200	1



BE/163

Beakers, tapered style, printed graduations, polypropylene

- Tapered stacking style - space saving, easy to handle
- High clarity PP - autoclavable at 121°C
- Durable and economical with crisp clear blue easy to read indelible graduations
- Precise non-dripping pouring spouts - safe and easy to use
- Full range of sizes, 25ml to 3000ml all with SPI recycling codes imprinted

Cat. No.	Code	Capacity ml	Sub-div. ml	Base x rim dia mm	Pk Qty
BE/16334	BPB0025P	25	5	32 x 41.5	10
BE/16340	BPB0050P	50	10	39 x 50.5	10
BE/16346	BPB0100P	100	20	48.5 x 61.5	10
BE/16356	BPB0250P	250	50	65.5 x 80.5	10
BE/16362	BPB0400P	400	50	74.5 x 92.5	10
BE/16366	BPB0500P	500	50	84 x 102.5	10
BE/16370	BPB0600P	600	50	84 x 104	5
BE/16376	BPB1000P	1000	100	99 x 120	5
BE/16380	BPB2000P	2000	200	125 x 163	1



BE/963

Beakers, polypropylene, three pouring spouts

Tricorn beakers with three pouring spouts and approximate moulded graduations.

Cat. No.	Description	Pk Qty
BE/96340	Beaker, PP, Tri-corn, 50ml	25
BE/96346	Beaker, PP, Tri-corn, 100ml	25
BE/96356	Beaker, PP, Tri-corn, 250ml	25
BE/96362	Beaker, PP, Tri-corn, 400ml	25
BE/96374	Beaker, PP, Tri-corn, 800ml	25
BE/96376	Beaker, PP, Tri-corn, 1000ml	25