

Bottles, reagent, narrow neck, glass stopper

- Pyrex borosilicate glass with ground neck
- Complies with ISO 4796-2

PYREX®

Cat. No.	Code	Capacity, ml	Neck	Pk Qty
BO/12702	1521/04D	100	14/15	10
BO/12704	1521/06D	250	19/26	10
BO/12706	1521/08D	500	24/29	10
BO/12708	1521/10D	1000	29/32	10
BO/12712	1521/16D	5000	45/40	1
BO/12714	1521/18D	10000	60/46	1
BO/12716	1521/22D	20000	60/46	1



BO/12708

Bottles, reagent, wide neck, glass stopper

Clear soda glass, Available in amber glass on request.

Cat. No.	Capacity, ml
BO/92202	50
BO/92204	100
BO/92206	250
BO/92208	500
BO/92210	1000



BO/922

Bottles, media lab, round shape, capped

Duran borosilicate glass bottles complete with cap and pouring ring complies with ISO 4796

**SCHOTT
DURAN®**

Cat. No	Code	Capacity, ml	Ht x Dia, mm	Thread, mm	Pk Qty
BO/22002	21801145*	25	70 x 36	25	10
BO/22004	21801175	50	88 x 46	32	10
BO/22006	21801245	100	100 x 56	45	10
BO/22008	21801365	250	138 x 70	45	10
BO/22010	21801445	500	176 x 86	45	10
BO/22012	21801545	1000	225 x 101	45	10
BO/22014	21801635	2000	260 x 136	45	10
BO/22016	21801735	5000	335 x 181	45	1
BO/22018	21801865	10000	410 x 227	45	1
BO/22020	21801885	15000	465 x 258	45	1
BO/22022	21801915	20000	510 x 300	45	1

*This small bottle features a glass lip and does not require a pouring ring
For caps and pouring rings see page 67



BO/22012

Bottles, media lab, round shape, amber, uncapped

Duran amber borosilicate glass bottles without cap or pouring ring

**SCHOTT
DURAN®**

Cat. No	Code	Capacity, ml	Ht x Dia, mm	Thread, mm	Pk Qty
BO/22222	2180614*	25	70 x 36	25	10
BO/22224	2180617	50	88 x 46	32	10
BO/22226	2180624	100	100 x 56	45	10
BO/22228	2180636	250	138 x 70	45	10
BO/22230	2180644	500	176 x 86	45	10
BO/22232	2180654	1000	225 x 101	45	10
BO/22234	2180663	2000	260 x 136	45	10
BO/22236	2180673	5000	335 x 181	45	1
BO/22238	2180686	10000	410 x 227	45	1
BO/22240	2180688	15000	465 x 258	45	1
BO/22242	2180691	20000	510 x 300	45	1

*This small bottle features a glass lip and does not require a pouring ring
For caps and pouring rings see page 67



BO/22232