



BO/462

Bottles, security, graduated, wide neck, PE

Will withstand considerable pressure before rupture. Used for transporting sensitive liquid and solid samples from location to location, secured by tagging points on bottle and polypropylene cap. Confidential notes can be placed in the space between insert sealing plug and cap.

Cat. No.	Capacity, ml	Grads., ml	Pk Qty
BO/46220	50	10	10
BO/46222	100	25	10
BO/46224	250	25	10
BO/46226	500	50	10
BO/46228	1000	100	5
BO/46230	2000	100	Each



BO/5630

Bottles, heavy duty, narrow neck, PE

Milk type bottle with polyethylene cap

Cat. No.	Capacity, ml	Pk Qty
BO/56302	250	10
BO/56304	500	10
BO/56306	1000	10



BO/5631

Bottles, with carrying handle, wide neck PE

Cat. No.	Capacity, litres
BO/56310	2.5
BO/56312	5
BO/56314	10



BO/5632

Bottles, with carrying handle, narrow neck PE

Cat. No.	Capacity, litres
BO/56320	2.5
BO/56322	5
BO/56324	10
BO/56326	25