

Bottles, plastic, LDPE

Bottles, plain round, LDPE

Natural bottles supplied without cap.

Cat. No.	Capacity, ml	Pk Qty
BO/96402	60	50
BO/96404	125	50
BO/96406	175	50
BO/96408	250	50

Caps for BO/964 series bottles, 22mm R6, Black, ribbed polypropylene with EPE wad.

Cat. No.	Description
BC/96222	Black cap with EPE liner, pack of 50

Bottles, round, narrow neck, LDPE

- Translucent, flexible leak proof LDPE bottles with PP screw caps
- Heavy duty and very durable in use.
- Ideal for a wide range of liquid storage, media preparation and sampling applications.
- Excellent general chemical resistance.
- Supplied in easy-to-handle shrink-wrapped packs.

Azlon

Cat. No.	Code	Capacity, ml	Pk Qty
BO/16202	BNL0030PN	30	10
BO/16204	BNL0060PN	60	10
BO/16206	BNL0100PN	100	10
BO/16208	BNL0125PN	125	10
BO/16210	BNL0150PN	150	10
BO/16212	BNL0250P	250	10
BO/16214	BNL0500P	500	10
BO/16216	BNL1000P	1000	5
BO/16218	BNL2000P	2000	5

Bottles, wide neck, LDPE

Translucent and flexible with excellent impact resistance. Ideal for shipping liquids. Leakproof.

NALGENE

Cat. No.	Code	Capacity, ml	Pk Qty
BO/36402	2103-0001	30	12
BO/36404	2103-0002	60	12
BO/36406	2103-0004	125	12
BO/36408	2103-0008	250	12
BO/36410	2103-0016	500	12
BO/36412	2103-0032	1000	6

Bottles, carboy, wide neck, LDPE

With heavy duty walls. For packaging and transporting. An additional external support is required for transporting corrosive materials. Also useful for storing and shipping chemicals and standards. Wide shoulder handles for easy carrying and pouring. Leakproof. Graduated in 5 litre increments.

NALGENE

Cat. No.	Code	Capacity, litre
BO/66310	2210-0020	10
BO/66315	2210-0040	15
BO/66320	2210-0050	20
BO/66325	2210-0065	25
BO/66350	2210-0130	50



BO/964



BO/162



BO/364



BO/663