



BW/564

Bottles, wash, LDPE

Bottles, wash, narrow neck, LDPE

General purpose wash bottle



Cat. No.	Capacity, ml	Pk Qty
BW/56430	250	5
BW/56432	500	5
BW/56434	1000	5



BW/1647

Bottles, wash, wide neck, plain, LDPE



- Square shoulder wash bottles featuring a wide neck for easy spill free filling.
- Fitted with colour coded one-piece PP closures that maintain their shape in use and allow easy identification.
- Closures feature precision jet-tip aperture for fine control.
- Bottles feature SPI recycling codes and volume information on their base.
- Supplied in shrink-wrapped packs.

Cat. No.	Code	Closure colour	Capacity, ml	Pk Qty
BW/16470	WGW230P	White	250	5
BW/16472	WGW231P	Red	250	5
BW/16474	WGW232P	Yellow	250	5
BW/16476	WGW233P	Green	250	5
BW/16478	WGW234P	Blue	250	5
BW/16480	WGW236P	White	500	5
BW/16482	WGW237P	Red	500	5
BW/16484	WGW238P	Yellow	500	5
BW/16486	WGW239P	Green	500	5
BW/16488	WGW235P	Blue	500	5
BW/16490	WGW443P	White	1000	5
BW/16492	WGW445P	Red	1000	5
BW/16494	WGW444P	Yellow	1000	5
BW/16496	WGW447P	Green	1000	5
BW/16498	WGW446P	Blue	1000	5



BW/1645

Bottles, wash, printed, narrow neck, LDPE



- Stylish labelled wash bottles enable easy identification of contents
- Indelibly printed with a choice of eight popular laboratory solvents, appropriate hazard symbols and CAS numbers
- The slim design enables comfortable handling and enhanced ease of use
- Fitted with colour coded PP, leak-proof, swivel dispensing closure
- Closures feature precision jet-tip aperture for fine control
- Supplied capped in shrink-wrapped packs of five

Cat. No.	Code	Printing and closure colour	Label	Capacity, ml
BW/16440	WGF530P	White	Water	250
BW/16442	WGF600P	White	Distilled Water	250
BW/16446	WGF531P	Red	Acetone	250
BW/16448	WGF532P	Green	Methanol	250
BW/16450	WGF533P	Blue	Isopropanol	250
BW/16452	WGF534P	Yellow	I.M.S.	250
BW/16454	WGF535P	Orange	Ethanol	250
BW/16456	WGF536P	White	Water	500
BW/16458	WGF601P	White	Distilled Water	500
BW/16462	WGF537P	Red	Acetone	500
BW/16464	WGF538P	Green	Methanol	500
BW/16466	WGF539P	Blue	Isopropanol	500
BW/16468	WGF540P	Yellow	I.M.S.	500
BW/16469	WGF541P	Orange	Ethanol	500