

Bottles, wash, printed, wide neck, LDPE

- Multilingual labelled wide neck wash bottles (English, French, German & Spanish)
- Indelibly printed with a choice of water types.
- Fitted with white one-piece PP closures that maintain their shape in use.
- Closures feature precision jet-tip aperture for fine control.

Azlon

Cat. No.	Code	Capacity, ml	Label	Pk Qty
BW/16424	WGW530PML	250	Water	5
BW/16426	WGW600PML	250	Distilled Water	5
BW/16428	WGW602PML	250	Deionised Water	5
BW/16324	WGW536PML	500	Water	5
BW/16328	WGW603PML	500	Deionised Water	5
BW/16330	WGW401PML	500	Sodium Hypochlorite	5



BW/163

Bottles, wash, venting, wide neck, LDPE

- Wide neck wash bottles featuring the patented Driplok vapour venting valve that prevents hazardous volatile solvent drips.
- Indelibly printed in 4 languages: English, French, German and Spanish.
- Choose from 5 solvent labels all including the Globally Harmonized Symbols, NFPA diamond, risk phrases and safety information.
- Fitted with colour coded one-piece PP closures that maintain their shape in use.
- Closures feature precision jet-tip aperture for fine control.
- Supplied capped in shrink-wrapped packs of five

Azlon

Cat. No.	Code	Label	Closure Colour	Pk Qty
250ml capacity				
BW/16350	WGW531VTML	Acetone	Red	5
BW/16352	WGW535VTML	Ethanol	Orange	5
BW/16354	WGW534VTML	I.M.S.	Yellow	5
BW/16356	WGW533VTML	Isopropanol	Blue	5
BW/16358	WGW532VTML	Methanol	Green	5
500ml capacity				
BW/16364	WGW537VTML	Acetone	Red	5
BW/16366	WGW541VTML	Ethanol	Orange	5
BW/16368	WGW540VTML	I.M.S.	Yellow	5
BW/16370	WGW539VTML	Isopropanol	Blue	5
BW/16372	WGW538VTML	Methanol	Green	5



BW/16358

Dispensing bottles, spouted cap, LDPE

- Round LDPE bottles fitted with spouted dropper closure.
- Closure features a captive sealing cap and can be cut to adjust flow rate.
- Suitable for use as a dropping or wash bottle.
- Good resistance to chemical attack.

Azlon

Cat. No.	Code	Capacity, ml	Pk Qty
BO/16430	BGE814PNC	30	10
BO/16432	BGE820PNC	60	10
BO/16434	BGE828PNC	150	10
BO/16436	BGE830PNC	250	10
BO/16438	BGE836PNC	500	10



BO/1643