

Bottles, graduated, square, wide neck, HDPE

Kartell LABWARE

Space saving bottles, ideal for long term storage of samples. Approximate graduations allow reproducible filling without the use of cylinders/beakers. Easy grip caps can be wired down for security. Head space between cap and insert plug can be used for confidential notes.

Cat. No.	Capacity, ml	Grads., ml	Pk Qty
BO/46382	25	5	10
BO/46384	50	10	10
BO/46386	100	20	10
BO/46388	250	50	10
BO/46390	500	100	10
BO/46392	1000	100	5
BO/46394	2000	100	Each



BO/463

Bottles, general purpose, wide neck, HDPE

NALGENE

Available in both natural and amber these general purpose bottles, with leakproof polypropylene screw closure, are durable and more rigid than LDPE. Suitable for shipping liquids.

Cat. No.	Code	Capacity, ml	Pk Qty
Natural			
BO/36124	2104-0001	30	12
BO/36126	2104-0002	60	12
BO/36128	2104-0004	125	12
BO/36130	2104-0008	250	12
BO/36132	2104-0016	500	12
BO/36134			
Amber			
	2104-0032	1000	6
BO/36140	2106-0002	60	12
BO/36142	2106-0004	125	12
BO/36144	2106-0008	250	12
BO/36146	2106-0016	500	12
BO/36148	2106-0032	1000	6



BO/3612

Bottles, round, wide neck, HDPE

Azlon

- Translucent, rigid and leak proof HDPE bottles with PP screw caps.
- Screw caps have non-stick feature and open and close in less than 1 turns. 1/2
- Heavy duty and very durable in use.
- Ideal for a wide range of liquid storage, media preparation and sampling applications.
- Excellent general chemical resistance.
- Supplied in shrink-wrapped packs.

Cat. No.	Code	Capacity, ml	Pk Qty
BO/16130	BWH0030PN	30	10
BO/16132	BWH0060PN	60	10
BO/16140	BWH0250PN	250	10
BO/16142	BWH0500PN	500	10
BO/16144	BWH1000PN	1000	5
BO/16146	BWH2000P	2000	5



BO/161

Bottles, square, wide neck, HDPE

NALGENE

Square shape saves space. With leakproof polypropylene screw cap.

Cat. No.	Code	Capacity, ml	Pk Qty
BO/36160	2114-0002	60	12
BO/36162	2114-0006	175	12
BO/36164	2114-0008	250	12
BO/36166	2114-0016	500	12
BO/36168	2114-0032	1000	6



BO/3616