



BO/4612

Bottles, wide neck, HDPE

Cylindrical bottles supplied with insert and cap. The very wide neck makes them especially suitable for powders and solids.

Cat. No.	Capacity, ml	Pk Qty
BO/46120	70	10
BO/46122	120	10
BO/46124	250	10
BO/46126	500	10
BO/46128	1000	5
BO/46130	2000	Each



BO/1612

Bottles, plastic, HDPE Bottles, square, graduated, wide neck



- Lightweight rectangular bottles with plug seal and screw cap.
- Moulded in graduations for approximate volume measurement.
- Ideal for powders and granules
- Space saving design.
- Supplied in shrink-wrapped packs.

Cat. No.	Capacity, ml	Pk Qty
BO/16120	500	10
BO/16122	1000	5
BO/16124	2000	5



BO/4610

Bottles, carboy, heavy walled, narrow neck, HDPE



Heavy walled (2mm) robust bottles which can withstand considerable abuse. The insert plug and soft seal cap makes them leakproof. The cap can be wired down for added safety. Complete with carrying handles.

Cat. No.	Code	Capacity, litres
BO/46102	1644	5
BO/46104	1646	10
BO/46106	1648	25
BO/46108	1650	50



BO/4611

Bottles, carboy, heavy walled, wide neck, HDPE

Robust bottle, ideal for liquids or powders. Wide neck for easy filling and emptying. Insert seal for completely leakproof closure. Caps can be wired down for safety. Tommy bar lugs on 10 litre capacity to assist release. Carrying handles.

Cat. No.	Capacity, litres
BO/46110	5
BO/46112	10