

Bottle, pistol grip, HDPE

Designed to allow comfortable grip and minimise fatigue. 500ml capacity bottle. Trigger spray adjustable from fine to single jet.

Cat. No.	Description	Pk Qty
BO/96110	Pistol grip bottle	10
BO/96112	Trigger spray	10

Bottle and trigger spray must be ordered separately.



BO/9611

Bottles, plastic, PP

Bottles, narrow neck, PP, Autoclavable

NALGENE®

Ideal for general lab applications. With leakproof polypropylene screw cap.

Cat. No.	Code	Capacity, ml	Pk Qty
BO/36202	2006-9125	4	12
BO/36204	2006-9025	8	12
BO/36206	2006-9050	15	12
BO/36208	2006-0001	30	12
BO/36210	2006-0002	60	12
BO/36212	2006-0004	125	12
BO/36214	2006-0008	250	12
BO/36216	2006-0016	500	12
BO/36218	2006-0032	1000	6



BO/362

Bottles, round, narrow neck, PP Autoclavable

Azlon

- Translucent, rigid and leak proof PP bottles with PP screw caps.
- Heavy duty and very durable in use.
- Ideal for a wide range of liquid storage, media preparation and sampling applications.
- Excellent general chemical resistance.
- Supplied in shrink-wrapped packs.

Cat. No.	Ref.	Capacity, ml	Height, mm	Pk Qty
BO/16220	BNP0030PN	30	67	10
BO/16222	BNP0060PN	60	79	10
BO/16224	BNP0100PN	100	89	10
BO/16226	BNP0125PN	125	103	10
BO/16228	BNP0150PN	150	117	10
BO/16230	BNP0250P	250	140	10
BO/16232	BNP0500P	500	160	10
BO/16234	BNP1000P	1000	200	5
BO/16236	BNP2000P	2000	237	5



BO/162