



BO/362

Bottles, wide neck, PP, Autoclavable

NALGENE®

For chemicals, specimens and general use. With leakproof polypropylene screw closure.

Cat. No.	Code	Capacity, ml	Pk Qty
BO/36220	2105-0001	30	12
BO/36222	2105-0002	60	12
BO/36224	2105-0004	125	12
BO/36226	2105-0008	250	12
BO/36228	2105-0016	500	12
BO/36230	2105-0032	1000	6



Natural

Bottles, round, wide neck, PP, Autoclavable

Azlon

- Translucent rigid and leak-proof PP bottles with PP screw caps.
- Available in natural and amber for light sensitive material
- Screw caps have nonstick feature and open or close in less than 1 turns for ease of ½ use.
- Heavy-duty and very durable in use.
- Autoclavable at 121°C.
- Ideal for a wide range of liquid or solids storage media preparation and sampling applications.
- Excellent general chemical resistance.
- Supplied in shrink-wrapped packs.

Cat. No.	Code	Capacity, ml	Height, mm	Pk Qty
Natural				
BO/16260	BWP0030PN	30	70	10
BO/16262	BWP0060PN	60	82	10
BO/16264	BWP0100PN	100	79	10
BO/16266	BLP0125P	125	100	10
BO/16268	BLP0150P	150	113	10
BO/16270	BWP0250P	250	141	10
BO/16272	BWP0500P	500	161	10
BO/16274	BWP1000N	1000	207	5
BO/16276	BWP2000P	2000	243	5
Amber				
BO/16240	BWP0060AP	60	84	10
BO/16242	BWP0150AP	150	104	10
BO/16244	BWP0250AP	250	141	10
BO/16246	BWP0500AP	500	161	10
BO/16248	BWP1000AP	1000	207	5
BO/16250	BWP2000AP	2000	243	5



Amber